



Oct 19, 2019

Western Blot (tank-blot) + Antibody staining

[iGEM Dusseldorf¹](#)¹Heinrich-Heine Universität Düsseldorf

1

Works for me

[dx.doi.org/10.17504/protocols.io.8gdhts6](https://doi.org/10.17504/protocols.io.8gdhts6)

iGEM Dusseldorf ⚡

MATERIALS TEXT

blotting buffer
SDS-Gel
100 % Ethanol
PDVF membrane

Fresh TBS-BSA
1x TBS-BSA, 100 µl Twin-20 and x µl Antibody
TBST

10x TBS:
24g Tris base
88 g NaCl
Dissolve in 900 ml H₂O
Adjust pH= 7.6 (HCl)
Fill to 1 L with H₂O
--> store at 4° C

1x TBS-BSA (prepare fresh):
50 ml 1x TBS
1.2 g BSA
--> store in fridge overnight

Tween-20 (50%):
10 ml Tween-20
10 ml water

TBST:
1 L 1x TBS
2 ml Tween-20
--> store at 4° C

1x TBS:
100 ml 10x TBS
900 ml H₂O
--> store at 4° C

Blotting buffer:
100 ml 10x SDS running buffer
200 ml 100% methanol
700 ml H₂O

SDS Running buffer (10x):
30 g TRIS
144 g Glycin
10 g SDS

Dissolve in 1 L H₂O

--> store at RT

Western blot

- 1 Create blotting buffer stock containing:
 - 100mL 10x SDS running buffer
 - 200mL 100% Methanol
 - 700mL H₂OFill a tray with blotting buffer
- 2 Transfer SDS gel to the blotting buffer
- 3 Assemble blotting sandwich (take it, soak sponge and place on it, soak filter paper and place on it, place inverted SDS gel on top without creating bubbles)
- 4 Add charged PDVF membrane with forceps (charged in 100 % ethanol and than soaked in buffer)
- 5 Add soaked Filter and Streak out bubbles
- 6 Add soaked sponge
- 7 Press together and close container
- 8 Place in blotting apparatus and add blotting buffer
- 9 Run for 1 h at 100 V 350-500 mA

Antibody staining

- 10 Add 40 ml fresh TBS-BSA in empty pipet tip box.
- 11 Place membrane from blotting apparatus with forceps into the solution
- 12 Protein side needs to be up
- 13 Incubate overnight at 4° C with shaking 20 rpm

- 14 Discard blocking solution
- 15 50 ml 1x TBS-BSA, 100 µl Tween-20 and x µl Antibody (dependend on the affinity)
- 16 Incubate for 1 h with shaking
- 17 Wash 2 times with TBST quickly , then wash 3 times a 10 min with TBST
- 18 If you are work with a secondary antibody icubate in the secondary antibody solution and repeat steps 17

Antibody staining

- 19 Develop using ECL (enhanced chemiluminescence) with HRP



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited