

1M Sodium citrate

Matthew Sullivan Lab

Abstract

Citation: Matthew Sullivan Lab 1M Sodium citrate. **protocols.io**

dx.doi.org/10.17504/protocols.io.c6azad

Published: 12 Oct 2015

Protocol

Step 1.

Add 294 g Na Citrate dihydrate (Mallinckrodt 0754-12) to 700ml MilliQ water

Step 2.

Stir to dissolve

Step 3.

Bring volume to 1L with MillQ water

Step 4.

Autoclave