

# Spore production for *Phytophthora infestans*

Remco Stam

## Abstract

**Citation:** Remco Stam Spore production for *Phytophthora infestans*. **protocols.io**

dx.doi.org/10.17504/protocols.io.fmmbk46

**Published:** 22 Aug 2016

## Protocol

### Step 1.

Take agar cubes / rains from glycerol stock and place 3-4 pieces on agar plates

When using another plate, cut small cubes (approx 0.5x0.5 cm) and place 3-4 upside down on the agar plate

For longer term cultivation use RyeA agar

For sporulation use RyeB agar

### PROTOCOL

#### . [Rye B Agar](#)

CONTACT: [Remco Stam](#)

#### Step 1.1.

Soak 60 g of rye grain in distilled water for 24 hours at room temperature. This is done in a small tray so that water just covers grain. Cover tray tightly with aluminum foil.

#### Step 1.2.

Next day, pour supernatant off germinated grain and put aside.

#### Step 1.3.

Boil rye grain (do not blenderize) for one hour in enough additional distilled water to cover the grain (about 1 inch above). Do this in a 2 liter beaker with foil to cover top. Check water level regularly, add more water if needed.

#### Step 1.4.

Filter through 4 thicknesses of cheese cloth, squeezing gently to remove residual liquid. Discard cheese cloth and grain

**Step 1.5.**

Combine the supernatants.

**Step 1.6.**

Add 15g of Agar, 20g of sucrose and 0.05g of  $\beta$ -Sitosterol. Adjust volume to 1l.

**Step 1.7.**

Autoclave at 15 psi for 20 min

**Step 2.**

Grow plates for about 7 days in the dark, until the plate is fully covered

Optimum temperature is about 18-20°C

**Step 3.**

Move plates to at 16/8 light/dark cycles for about 2-3 days to induce sporulation

**Step 4.**

Harvest sporangia by adding several ml ice cold water to the plates and scraping of the sporangia with a pipette tip

**Step 5.**

Aspirate all liquid and store in fridge for 3-6 hours to induce hatching of the spores

**Step 6.**

Pipette off the upper layers (spores swim up, debris sink) and check spore density under microscope