

# E. coli Heat Shock

Alex Rajewski

## Abstract

**Citation:** Alex Rajewski E. coli Heat Shock. **protocols.io**

dx.doi.org/10.17504/protocols.io.nwfdfbn

**Published:** 19 Mar 2018

## Protocol

Set up

### Step 1.

Thaw competent cells on ice.

Set up

### Step 2.

Warm SOC media to room temperature.

Set up

### Step 3.

Heat a dry block to 42°C

Set up

### Step 4.

Label a 1.5mL eppendorf tube for each transformation reaction with the construct name and place it on ice.

Heat Shock

### Step 5.

Add 2µL of plasmid to the labeled eppendorf tube on ice.

 AMOUNT

2 µl Additional info: plasmid

Heat Shock

### Step 6.

Add 45µL of thawed competent cells to the tube with the plasmid and gently swirl with the pipet tip to mix.

## AMOUNT

45 µl Additional info: competent cells

### Heat Shock

#### Step 7.

Incubate the mixture on ice for 30 minutes.

### Heat Shock

#### Step 8.

Place the mixture in the dry block at 42°C for 30 seconds

### Heat Shock

#### Step 9.

Incubate the mixture on ice for 2 minutes.

### Grow Up

#### Step 10.

Add 250µL of SOC to the mixture and shake for 1 hour at 37°C and 250rpm

### Grow Up

#### Step 11.

Spin for 30 seconds at 6000rpm to pellet

## GOTO

Concentrate -> go to step #1

### Grow Up

#### Step 12.