

# AMPure purification protocol

Tomasz Suchan

## Abstract

**Citation:** Tomasz Suchan AMPure purification protocol. **protocols.io**

dx.doi.org/10.17504/protocols.io.re5d3g6

**Published:** 16 Jul 2018

## Guidelines

The protocol below works for 200 µl PCR tubes or plates. If working with 1.5 ml tubes, use larger volumes of EtOH for the washes (step 6 and 7) and Tris/water for the final elution (step 9).

## Protocol

### Step 1.

Remove beads from the fridge, equilibrate to the room temperature, mix thoroughly but do not vortex.

### Step 2.

Add desired ratio of AMPure beads to the purified sample and mix well by pipetting.

### Step 3.

Incubate 5 minutes.

### Step 4.

Place on the magnetic rack.

### Step 5.

Let it stand for 5 minutes on the rack, aspirate and discard supernatant.

### Step 6.

Add 200 µl of freshly prepared 70% EtOH, incubate 30 seconds, aspirate and discard EtOH.

### Step 7.

Repeat the wash: add 200 µl of freshly prepared 70% EtOH, incubate 30 seconds, aspirate and discard EtOH.

### Step 8.

Wait until the pellet dries out completely. DO NOT remove the plate/tubes from the magnetic rack, as the dried AMPure might pop out from the tubes!

### Step 9.

Add desired volume of Tris 10 mM or water, wetting the dried AMPure pellets (add 1 µl to the final volume to avoid pipetting out the beads)

---

**Step 10.**

Remove from the magnetic rack.

---

**Step 11.**

Resuspend by pipetting or vortexing.

---

**Step 12.**

Incubate 10 minutes, incubating in 37°C can improve DNA yield.

---

**Step 13.**

Place on the magnetic rack.

---

**Step 14.**

Let it stand for 5 minutes, pipette out and save supernatant. The eluted DNA is in the supernatant, do not discard it!

---