Restriction Digest Reaction Version 2

New England Biolabs

Abstract

Citation: New England Biolabs Restriction Digest Reaction. protocols.io

dx.doi.org/10.17504/protocols.io.iszcef6

Published: 02 Mar 2018

Protocol

Step 1.

DNA 1µg



1 μg Additional info:

Step 2.

10X NEBuffer **5μl** (1X)

AMOUNT

5 μl Additional info:

Step 3.

Restriction Enzyme, 10 units is sufficient, generally 1µl is used