



Jul 12, 2019

Mouse Stellate Immunohistochemistry protocol

Pradeep Rajendran¹

¹[University of California, Los Angeles]

1

Works for me

dx.doi.org/10.17504/protocols.io.2qggdtw



scott john



- 1 See mouse stellate isolation protocol
- 2 Complete/intact stellate ganglia are fixed in 4% paraformaldehyde overnight at 4°C.
- 3 Fixed tissue is rinsed in phosphate buffered saline (PBS), and stored in PBS + 0.02% sodium azide.
- 4 The stellate ganglion whole mounts are permeabilized in 'block' solution. (10% normal donkey serum in 1.0% BSA + 0.4% Triton X-100 + PBS).
- 5 Ganglia are stained with antibodies directed to either tyrosine hydroxylase (TH) or protein gene product 9.5 (PGP 9.5). Secondary staining with streptavidin conjugated to ATTO-647N was used to visualize neurobiotin filling. Stained tissue was rinsed in PBS and mounted on glass slides.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited