

# Fixation and dehydration protocol for Scanning Electron Microscopy (SEM), for the observation of morphology of secretory spines

Mario A. Sandoval-Molina, Hilda Araceli Zavaleta-Mancera, Simón Morales-Rodríguez, Mariusz Krzysztof Janczur

## Abstract

The presence and structure of EFNs in *Opuntia robusta* had not been investigated. We used light, scanning-electron, and transmission-electron microscopy to examine morphology, anatomy, and ultrastructure of the secretory spines in areoles in female and hermaphrodite individuals of *O. robusta*. Young cladodes develop areoles with modified and secretory spines as EFNs only active during the early growth phase of female and hermaphrodite individuals. EFNs are non-vascularized structures, with no stomata, that consist of three distinct tissues: a basal meristematic tissue; a middle elongation region; and an apical secretory cone (asc) formed by large globular epidermal cells, of sac shape, containing nectar and medullar elongated cells.

**Citation:** Mario A. Sandoval-Molina, Hilda Araceli Zavaleta-Mancera, Simón Morales-Rodríguez, Mariusz Krzysztof Janczur Fixation and dehydration protocol for Scanning Electron Microscopy (SEM), for the observation of morphology of secretory spines. **protocols.io**

dx.doi.org/10.17504/protocols.io.k3vcyn6

**Published:** 01 Dec 2017

## Protocol

### Stock solutions

#### Step 1.

Prepared with deionised water.

A) 0.2 M of  $\text{NaH}_2\text{PO}_4$  (Sodium phosphate monobasic anhydrous, CAS Number: [7558-80-7](#)) BioUltra,  $\geq 99.0\%$  (T) (Sigma-Aldrich S3139).

B) 0.2 M of  $\text{Na}_2\text{HPO}_4$  (Sodium phosphate dibasic anhydrous, CAS Number: 7558-79-4) Ultra  $\geq 99.5\%$ , (Sigma-Aldrich S7907).

To make 200 mL of phosphate buffer solution (0.1 M, pH 7.2), mix 28 mL of stock solution A + 72 mL of stock solution B and 100 mL of deionised water. Adjust pH with solutions of 1N HCl or 5 M NaOH if necessary.

### Buffered glutaraldehyde solution

## **Step 2.**

### **To prepare 100 mL**

Mix 10 mL of glutaraldehyde solution (25% for electron Microscopy, Merck, Germany) with 90 ml of phosphate buffer Sorensen's solution (0.1 M, pH 7.2).

### **1% Osmium tetroxide in deionised water**

Dissolve 0.25 g of  $\text{OsO}_4$  in 25 mL of deionised water (use fume hood). Protect the solution from light and store at 4 °C.

## **Scanning Electron Microscopy (SEM)**

### **Step 3.**

### **Critical/point dried and SEM observations.**

Samples were critically-point dried using a Samdri-780 (Tousim Research Corporation, Rockville, USA), mounted on cooper studs with a carbon conductive tape, double coated (Ted Pella I.N.C U.S.A) and then sputter-coated with gold-palladium with a Fine coat ion sputter (JFC-1100, Jeol Ltd., Tokyo, Japan). Observations were performed using a JEOL JSM 6390 scanning electron microscope operated at 25 kV.