

Quantification of circulating microRNA using single TaqMan assays

Helle Glud Binderup, Jonna Skov Madsen, Kim Houllind, Rikke Fredslund Andersen, Claus Lohman Brasen

Abstract

A protocol for quantification of circulating microRNA using TaqMan assays. MicroRNAs are purified from platelet-poor plasma using Nucleospin columns. During the purification step samples are spiked with cel-miR-39 as a mean of technical normalization. Normalization is performed using the ΔC_t -method.

Citation: Helle Glud Binderup, Jonna Skov Madsen, Kim Houllind, Rikke Fredslund Andersen, Claus Lohman Brasen. Quantification of circulating microRNA using single TaqMan assays. **protocols.io** dx.doi.org/10.17504/protocols.io.q9edz3e

Published: 05 Jul 2018

Guidelines

Blood samples should be obtained using a minimum of venous stasis and with discard of the first 3 mL of blood.

Platelet-poor plasma should be prepared within 2 hours from blood sampling

Before start

We recommend to use a few samples to test if the volume of spike-in added during microRNA purification is adequate before purifying all your samples.

The C_t -values of the spike-in and the target miRNA should be within the same range, otherwise adjust the volume of spike-in added during miRNA purification.

Materials

K2-EDTA containing tubes 366643 by Becton-Dickinson

Nucleospin® miRNA Plasma 740971.50 by Macherey-nagel

TaqMan® MicroRNA Reverse Transcription Kit 4366597 by Applied Biosystems, Foster City, California, USA

2X TaqMan Universal PCR Master Mix

4318157 by Applied Biosystems, Foster City,
California, USA

20X TaqMan MicroRNA Assay 4440887 by
Thermo Fisher Scientific

Protocol

Preparation of platelet-poor plasma (PPP)

Step 1.

Dual centrifugation

1. use 10 ml of EDTA anticoagulated whole blood
2. centrifugation at 3000 g for 15 minutes (acceleration 5, brake 6, temperature 18 °C)
3. transfer plasma phase to new tube, leaving approximately 1 mL of plasma on top of the buffy coat
4. centrifugation at 3000 g for 15 minutes (acceleration 5, brake 6, temperature 18 °C)
5. transfer plasma phase to cryo tubes, leaving approximately 1 ml of plasma in the bottom of the tube
6. store at -80 °C

NOTES

Helle Glud Binderup 22 Jun 2018

Alternatively, a prolonged single centrifugation may be used

1. use 5 ml of EDTA anticoagulated whole blood
2. centrifuge at 3000 g for 30 minutes (acceleration 5, brake 6, temperature 18 °C)
3. transfer plasma phase to cryo-tubes, leaving approximately 0.5 ml of plasma on top of the buffy coat

EQUIPMENT

Equipment brand:
Hettich centrifuge
SKU:
4706-01
Specifications:
Rotina 420R

MicroRNA purification















Step 2.

REAGENTS

Nucleospin® miRNA Plasma
740971.50 by Macherey-nagel

Step 3.

Follow the instructions given by the manufacture (se notes before you start):

NucleoSpin® miRNA Plasma		
1 Prepare sample		300 µL plasma or serum* 90 µL MLP Vortex 5 s RT, 3 min
2 Precipitate protein	 	30 µL MPP Vortex 5 s RT, 1 min 11,000 x g, 3 min
3 Transfer supernatant		Transfer clear supernatant to Collection Tube (2 mL, lid)
4 Adjust binding conditions	 	400 µL isopropanol Vortex 5 s
5 Bind RNA and DNA	 	Load sample on NucleoSpin® miRNA Column RT, 2 min 11,000 x g, 30 s
6 Optional: Digest DNA	 	Optional: 1st 700 µL MW2 11,000 x g, 30 s 2nd 250 µL MW2 11,000 x g, 2 min 50 µL rDNase in Reaction Buffer for rDNase RT, 15 min
7 Wash and dry	 	1st 100 µL MW1 11,000 x g, 30 s 2nd 700 µL MW2 11,000 x g, 30 s 3rd 250 µL MW2 11,000 x g, 2 min
8 Elute RNA	 	30 µL RNase-free H ₂ O RT, 1 min 11,000 x g, 1 min

📌 NOTES

Helle Glud Binderup 22 Jun 2018

STEP 1: use 300 µL platelet-poor plasma

STEP 4: after addition of isopropanol, add also 5 µL spike-in (cel-miR-39, 2.75×10^{-12} M)

STEP 6: perform the optional DNA digest

Reverse transcription

Step 4.

REAGENTS

TaqMan® MicroRNA Reverse Transcription Kit 4366597 by Applied Biosystems, Foster City, California, USA

20X TaqMan MicroRNA Assay 4440887 by Thermo Fisher Scientific

Reverse transcription

Step 5.

Prepare RT master mix:

Component	Master mix volume per 15 µL reaction*
100 mM dNTPs (with dTTP)	0.15 µL
Multiscribe RT enzyme (50 U/µL)	1 µL
10x RT buffer	1.5 µL
RNase inhibitor (20 U/µL)	0.19 µL
Nuclease free water	to 15 µL in total**

*add 10-20% excess volume

**add water to a total reaction volume of 15 µL (including microRNA and RT-primers, see step 7)

Mix gently and place on ice

NOTES

Helle Glud Binderup 25 Jun 2018

Remember to include a RT-negative sample (no template)

It is also a good idea to include a RT-positive sample (a microRNA-sample included in all runs)

Step 6.

For each RT reaction, combine RT master mix with 2 µL of purified microRNA in a 0.2 µL polypropylene reaction tube

(the volume of RT master mix is dependent on the number of RT-primers, see step 7)

Keep on ice

Step 7.

Prepare RT primer mix:

Add for each RT reaction 0.75 µL 20x RT primer from each microRNA assay set

Example with two microRNA assays:

Component	Volume per 15 µL reaction*
20x RT specific primer #1	0.75 µL
20x RT specific primer #2	0.75 µL
Total volume	1.5 µL**

*add 10-20% excess volume

**when using two microRNA assays the volume of water needed in step 5 is 8.66 µL for each reaction, and the volume of master mix added in step 6 is 11.5 µL

Mix gently

Step 8.

Add the appropriate volume of RT primer mix to the reaction tubes prepared in step 6

Mix gently and incubate on ice for 5 minutes

Step 9.

Transfer the reaction tubes to a thermocycler

Incubate at 16 °C for 30 min, at 42 °C for 30 min and at 85 °C for 5 min.

Step 10.

Continue immediately to the PCR amplification or store the RT-reaction in -20 °C.

Quantitative real-time PCR

Step 11.



REAGENTS

2X TaqMan Universal PCR
Master Mix 4318157 by Applied
Biosystems, Foster City, California,
USA

20X TaqMan MicroRNA Assay
4440887 by Thermo Fisher Scientific

Quantitative real-time PCR

Step 12.

Prepare qPCR reaction mix (run each sample in doublets):

Component	Volume per 20.3 µL reaction*
2x TaqMan Universal PCR master mix	10 µL
20x TaqMan microRNA assay	1 µL
Nuclease free water	8 µL
Total volume	19 µL

*add 10-20% excess volume

Mix gently

Quantitative real-time PCR

Step 13.

Combine 19 µL of qPCR reaction mix with 1.3 µL of RT-reaction (cDNA) in a 96 well plate
Seal plate

Quantitative real-time PCR

Step 14.

 SOFTWARE PACKAGE

Sequence Detection System (SDS), 2.4

<https://www.thermofisher.com/search/results?query=sequence+analysis+software&persona=Catalog&navId=10949&refinementAction=true&focusarea=Search+All>

Quantitative real-time PCR

Step 15.

Transfer plate to ABI Prism 7900HT

Use the SDS software to set up the run

Incubate in 50 °C for 2 min and 95 °C for 10 min

Proceed with 40 cycles of 95 °C for 15 sec and 60 °C for 60 sec.

Normalization

Step 16.

Normalization is performed using the ΔC_t -method ($2^{-\Delta C_t}$)

$\Delta C_t = \text{mean } C_{t_{\text{target miRNA}}} - \text{mean } C_{t_{\text{cel-miR-39}}}$ (mean of doublets)

