

Mar 01,  
2019

Working

## *E. coli* K12 DNA Extraction

Version 3

Kenneth Schackart<sup>1</sup>, Kattika Kaarj<sup>1</sup>

<sup>1</sup>University of Arizona

dx.doi.org/10.17504/protocols.io.yrpfv5n

 Kenneth Schackart 

### ABSTRACT

How to extract DNA from *E. coli* K12 using Wizard® Genomic DNA Purification Kit by Promega®.

I do not claim any credit for the development of this protocol. It has been adapted from the protocol detailed in:

 [Wizard Genomic DNA Purification.pdf](#)

### PROTOCOL STATUS

#### Working

We use this protocol in our group and it is working

### MATERIALS



NAME	CATALOG #	VENDOR
Wizard(R) Genomic DNA Purification Kit	A1620	Promega

### MATERIALS TEXT



Additional materials:

- 1.5 mL microcentrifuge tubes
- Isopropanol, room temperature
- 70% ethanol, room temperature
- 




#### Culture bacteria

- 1 Culture *E. coli* K12 in BHI broth overnight.  
 2 mg lyophilized *E. coli* K12 in  10 ml BHI broth.





#### Pellet the cells

- 2 Add  1 ml cell suspension to 1.5 mL microcentrifuge tube.
- 3 Centrifuge at 13,000-16,000 × *g* for  00:02:00 .
- 4 Remove supernatant.





#### Lyse nuclei

- 5 Add  600  $\mu$ l of Nuclei Lysis Solution.
- 6 Gently pipet until the cells are resuspended.
- 7 Incubate at  80 °C on heating block for  00:05:00 to lyse the cells.
- 8 Cool to room temperature.

#### Degrade RNA

- 9 Add  3  $\mu$ l RNase Solution to the cell lysate.
- 10 Invert 2-5 times to mix.
- 11 Incubate at  37 °C for  00:15:00 to  01:00:00 .
- 12 Cool to room temperature.

#### Precipitate proteins

- 13 Add  200  $\mu$ l of Protein Precipitation Solution to the RNase-treated cell lysate.
- 14 Vortex vigorously at high speed for  00:00:20 .
- 15 Incubate on ice for  00:05:00
- 16 Centrifuge at 13,000-16,000  $\times g$  for  00:03:00 .

#### Harvest DNA




- 17 Transfer the supernatant containing the DNA to a clean 1.5 mL microcentrifuge tube containing  600  $\mu$ l isopropanol.







Some supernatant may remain in the original tube containing the protein pellet. Leave this residual to avoid contaminating the DNA solution with the precipitated protein.

- 18 Gently mix by inversion until the thread-like strands of DNA form a visible mass.

#### Wash and dry DNA

- 19 Centrifuge at 13,000-16,000 × *g* for  00:02:00 .
- 20 Carefully pour off the supernatant and drain the tube on clean absorbent paper.
- 21 Add  600 µl of 70% ethanol and gently invert the tube several times to wash the DNA pellet.
- 22 Centrifuge at 13,000-16,000 × *g* for  00:02:00 .
- 23 Carefully aspirate the ethanol.
- 24 Drain the tube on clean absorbent paper and allow to air-dry for 10-15 minutes.

#### Rehydrate DNA

- 25 Add  100 µl of DNA rehydration solution to the tube.
- 26 Rehydrate by incubating the solution overnight at room temperature or  4 °C .
- 27 Store DNA at  2 °C to  8 °C .



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited