

Human subchondral osteoblasts cell culture

Christelle Sanchez

Abstract

Method of isolation and culture of osteoblast coming from human subchondral bone of the knee

Citation: Christelle Sanchez Human subchondral osteoblasts cell culture. **protocols.io**

dx.doi.org/10.17504/protocols.io.mkmc4u6

Published: 12 Jan 2018

Materials

HEPES H6147 by [Sigma Aldrich](#)

Trypsin EDTA 25-051-Cl. by [Gibco - Thermo Fischer](#)

200 mM L-Glutamine G7513 by [Sigma](#)

DMEM 4.5 g/L glucose BE12-614F by [Lonza](#)

Collagenase IA C9891 by [Sigma Aldrich](#)

✓ Fetal Bovine Serum S1810 by Contributed by users

Penicilin Streptomycin DE17-602E by [Lonza](#)

PBS 1x BE17-516F by [Lonza](#)

Protocol

Subchondral bone dissection

Step 1.

Tibial subchondral bone plates were collected.

After careful elimination of trabecular bone and articular cartilage, subchondral bone was dissected to separate sclerotic from non-sclerotic zones. Non-sclerotic and sclerotic bone zones were identified by a marked difference in thickness. We considered sclerotic bone to be only that from the subchondral areas of bone with a thickness >2 mm and being either denuded or overlaid by fibrillated cartilage. We considered nonsclerotic bone to be only that from the subchondral areas of bone with a maximal thickness of 1 mm. Intermediate zones of the subchondral bone plate were discarded.

The plates were then cut into small fragments of approximately 2 mm³ with a bone rongeur, washed with DMEM and then submitted to enzymatic digestions.

Subchondral bone digestion

Step 2.

Small pieces of bone (2 mm³) were sequentially incubated with 1 mg/ml clostridial collagenase IA (Sigma-Aldrich) for 35 and 240 min successively (2 g of bone in 20 ml of enzymatic solution).

Subchondral osteoblast culture

Step 3.

The digested bone pieces were extensively rinsed in DMEM, placed into T-75 flasks (0.5 to 1 g of bone per flask for not sclerotic and 1.5 to 2 g of bone for sclerotic) and cultured in DMEM supplemented with 15% fetal bovine serum (FBS), 10 mM HEPES, 100 U/ml penicillin and 100 µg/ml streptomycin, until osteoblasts migrated out of bone explants (one week). Media was replaced twice a week. At this point, the medium was replaced with fresh media containing 10% FBS, 10 mM HEPES, 100 U/ml penicillin, 100 µg/ml streptomycin and 2 mM glutamine. At confluence, cells were collected by trypsinization, seeded (20,000 cells/cm²) in 6, 12 or 24-well plates (NUNC Nunclon plates) and grown for three days in DMEM containing 10% FBS, 100 U/ml penicillin, 100 µg/ml streptomycin, 10 mM HEPES and 2 mM glutamine.

Subchondral osteoblast culture

Step 4.