



Jun 26, 2019 Cellular protein extraction and Western blotting using dry transfer (iBlot system)

Ralitsa R. Madsen<sup>1</sup>

<sup>1</sup>University of Edinburgh

Working dx.doi.org/10.17504/protocols.io.4r4gv8w



## ABSTRACT

A generic protein extraction and Western blotting protocol, allowing both for total cell lysis or cytoplasmic protein extraction only. The cytoplasmic extraction version has been used when Western blotting for PI3K and MAPK/ERK signalling components; total cell lysis with sonication has been used when Western blotting for histones.

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

This protocol is an amalgamate of previous protocols provided by Dr Gemma Brierley (Institute of Metabolic Science, University of Cambridge) and Dr Tijana Mitic (CVS, University of Edinburgh), with additional modifications based on my own PhD/postdoc experience.

#### **GUIDELINES**

The protocol assumes that the cells have been collected as follows prior to lysis:

- 1) Washed once with cold DPBS.
- 2) Snap-frozen on dry ice or liquid nitrogen, followed by long-term storage at  $\,$  8 -80  $^{\circ}\text{C}$  .

#### MATERIALS

NAME Y	CATALOG #	VENDOR V
cOmplete™, EDTA-free Protease Inhibitor Cocktail	05056489001	Sigma Aldrich
20X MES Buffer	NP0002	Thermo Fisher Scientific
NuPAGE Antioxidant	NP0005	Thermo Fisher Scientific
Molecular Biology Grade Water	10154604	Fisher Scientific
BSA	A7906	Sigma Aldrich
NUPAGE LDS sample buffer (4x)	NP0007	Thermo Fisher Scientific
1x NUPAGE MOPS SDS running buffer (20x)	NP0001	Thermo Fisher Scientific
DC™ Protein Assay Kit	500-0112	BIO-RAD
cOmplete ULTRA Tablets Mini EasyPack PhosStop	4906845001	Sigma Aldrich
Precision Plus Protein™ Dual Color Standards	1610374	BIO-RAD
NuPAGE Sample Reducing Agent (10X)	NP0009	Thermo Fisher Scientific
NuPAGE 4-12%BT midi 12 2 well PAGE gels	WG1401BOX	Thermo Fisher Scientific
iBlot™ 2 Transfer Stacks nitrocellulose regular size	IB23001	Thermo Fisher Scientific
Immobilon Western Chemiluminescent HRP Substrate	WBKLS0500	Sigma Aldrich
EZBlue™ Gel Staining Reagent	G1041	Sigma Aldrich
Fisherbrand™ Cell Scrapers	11587692	Fisher Scientific
XCell4 SureLock™ Midi-Cell	WR0100	Thermo Fisher Scientific

MATERIALS TEXT

TBS/T for washes; [M]20 Milimolar (mM) Trizma base, [M]150 Milimolar (mM) NaCl, pH = 7.6 (tip: prepare a larger volume of 10X stock solution and dilute each time; use deionised H<sub>2</sub>O).

Some reagents are toxic if inhaled or similar and should be handled according to the accompanying materials safety sheet.

```
BEFORE STARTING
```

### Prepare the required lysis buffers depending on application

The buffers can be prepared in bulk and stored at 8.4 °C (we have previously stored buffers A and B for 2-3 years, amd they continued to work). Always check for precipitates in the stock before preparing a working aliquot.

Take out smaller aliquots when required and supplement with the indicated inhibitors. It might be a good idea to sterile-filter the working aliquot before use.

\* Add fresh before use.

Type A buffer: whole-cell lysis; relatively mild (used for immunoprecipitation of phosphoproteincs and regular WB of cytoplasmic proteins)

```
[M120 Milimolar (mM) HEPES (pH = 7.4)
[M1150 Milimolar (mM) NaCl
[M11.5 Milimolar (mM) MgCl2
[M10 Volume Percent Glycerol
[M11 Volume Percent TritonX-100
[M11 Milimolar (mM) EGTA (dissolves at & 37 °C and pH 8, stirring required)
```

Store at § 4 °C - when needed, prepare a 10 ml working aliquot and supplement with the following immediately before use:

[M]1 Milimolar (mM) \* PMSF (prepared as 100X stock dissolved in isopropanol and store as smaller aliquots at 8 -20 °C; short half-life once in aqueous solution)

[M]2 Milimolar (mM) \* Na<sub>3</sub>VO<sub>4</sub> (has to be dissolved in a particular way, for more details follow this link; prepare as 100X stock solution and store as smaller aliquots at 8 -20 °C)

- \* 1X EDTA-free protease inhibitor tablet (to 10 ml buffer)
- \* 1X PhosStop tablet (to 10 ml buffer)

Type B buffer: whole-cell lysis; relatively mild (used for nitrocellulose-based reverse phase protein arrays and regular WB of cytoplasmic proteins; for original reference, see Macleod et al. 2017 doi:10.1007/978-1-4939-7201-2)

```
[M]150 Milimolar (mM) HEPES (pH = 7.4)

[M]150 Milimolar (mM) NaCl

[M]1.5 Milimolar (mM) MgCl<sub>2</sub>

[M]10 Volume Percent Glycerol

[M]1 Volume Percent TritonX-100

[M]1 Milimolar (mM) EGTA (dissolves at & 37 °C and pH 8, stirring required)

[M]100 Milimolar (mM) NaF

[M]10 Milimolar (mM) Na<sub>4</sub>P<sub>2</sub>O<sub>7</sub>
```

Adjust pH to 7.4 and store at 8.4 °C - when needed, prepare a 10 ml working aliquot and supplement with the following immediately before use:

[M]2 Milimolar (mM) \* Na<sub>3</sub>VO<sub>4</sub> (has to be dissolved in a particular way, for more details follow this link; prepare as 100X stock solution and store as smaller aliquots at & -20 °C)

- \* 1X EDTA-free protease inhibitor tablet (to 10 ml buffer)
- \* 1X PhosStop tablet (to 10 ml buffer)

Type C buffer: whole-cell lysis; harsh for nuclear lysis (used for WB of histones) - modified RIPA buffer (higher SDS concentration used)

```
[M]20 Milimolar (mM) Tris-HCl (pH = 7.5)
[M]150 Milimolar (mM) NaCl
[M]1 Milimolar (mM) Na<sub>2</sub>EDTA
```

```
[M] 1 Milimolar (mM) EGTA (dissolves at § 37 °C and pH 8, stirring required)
   [M]1 Mass/Volume Percent SDS
   [M]1 Mass/Volume Percent Sodium deoxycholate
   [M]2.5 Milimolar (mM) Na<sub>4</sub>P<sub>2</sub>O<sub>7</sub>
   [M] 1 \text{ Milimolar (mM)} \beta-glycerophosphate
  1 μg/ml leupeptin
 Store at 🐧 4 °C - when needed, prepare a 10 ml working aliquot and supplement with the following immediately before use:
   [M]1 Milimolar (mM) * Na<sub>3</sub>VO<sub>4</sub> (has to be dissolved in a particular way, for more details follow this link; prepare as 100X stock
 solution and store as smaller aliquots at 8 -20 °C)
 * 1X EDTA-free protease inhibitor tablet (to 10 ml buffer)
 * 1X PhosStop tablet (to 10 ml buffer)
Cell lysis and protein extraction
      Pre-chill the benchtop microcentrifuge to § 4 °C .
      Prepare cell scraper and a beaker with PBS - this will be used to rinse the cell scraper in between processing of different samples.
 2
      Aliquot 10 ml of stock lysis buffer and dissolve the required supplements that need to be added fresh (see "Guidelines" for buffer details).
       Keep on ice.
 3
      Allow the snap-frozen cells to thaw on ice (one plate at a time, process each well for scraping and rinse the scraper in PBS in between
 4
      Add \[ \] 150 \muI of ice-cold lysis buffer per well and scrape the cells on ice.
               The volume can be reduced to 100 µl depending on expected yield. I use this volume when working with human pluripotent
               stem cells which give a high protein yield from a 6-well (3-5 mg/ml); this may not be the case with other cell lines and should
               be tested empirically.
      Transfer cell lysates to pre-labelled and pre-chilled tubes.
                                                                                                                                                                              30m
 6
      Incubate for minimum (§ 00:30:00 on ice.
                                                                                                                                                                              05s
 7
      Vortex each sample for © 00:00:05.
                                                                                                                                                                               5m
 8
      *
      If performing nuclear lysis for extraction of histones: sonicate the samples on a Diagenode Bioruptor using 5 pulses of 3 © 00:00:30 ON,
       © 00:00:30 OFF; setting =high. Following sonication, the lysates should appear clear/runny.
      Centrifuge the lysate at § 4 °C and (3)12000 x g for (5)00:10:00.
```

- 9.1 In the mean time, thaw previously prepared BSA standards (0 2 mg/ml).
- 10 Transfer the supernatant to new pre-chilled tubes and discard the pellets.



Stopping point: the samples can be stored at 8-20 °C at this point and processed for protein concentration measurements at a later time. Freeze-thawing should generally be limited, but up to 3 times has worked well for me in the past (alternatively, make multiple aliquots to avoid freeze-thawing of the entire volume each time).

11 Dilute the samples for protein concentration measurements (usually 1:5 but will depend on the exact samples and the expected yield). Use molecular-grade H<sub>2</sub>O.



Remember to dilute the lysis buffer the same.

It is easiest to perform these and all subsequent dilutions in PCR strip tubes.

- 12 Prepare DC assay solution A+S according to the manufacturer's instructions (BioRad): use 20 μl Reagent S to 11 ml Reagent A (prepare a mastermix in excess of what you will need; e.g. for loading of 8 standards and 8 samples in quadruplicate = 25 \* 16 \* 4 \* 1.2 (excess) = 1920 ~ 2 ml (2 ml Reagent A + 40 μl Reagent S).
- 13 Add 25 μl of the A+S mix to each well of a 96-well plate for protein concentration measurements (e.g. Sterilin™ Clear Microtiter™ Plates; Thermo Scientific Sterilin 611F96 (for absorbance measurements) #11349163)
- 14 Add 35 μl of each sample or standard dilution in triplicate; use multichannel pipettor if dealing with multiple samples and all diluted in PCR strip tubes.

Example layout

Plate 1	1	2	3	4	5	6	7	8	9	10	11	12
A	S1 (0 mg/ml)	S1	S1	G7_25	G7_25	G7_25	Lysis buffer dilution	Lysis buffer dilution	Lysis buffer dilution			
В	S2 (0.125 mg/ml)	S2	S2	G7_26	G7_26	G7_26						
С	S3 (0.250 mg/ml)	S3	S3	G7_27	G7_27	G7_27						
D	S4 (0.5 mg/ml)	S4	S4	G7_28	G7_28	G7_28						
E	S5 (0.750 mg/ml)	S5	S5	G7_29	G7_29	G7_29						
F	S6 (1 mg/ml)	S6	S6	G7_30	G7_30	G7_30						
G	S7 (1.5 mg/ml)	S7	S7	G7_31	G7_31	G7_31						
Н	S8 (2 mg/ml)	S8	S8	G7_32	G7_32	G7_32						

I prepare my BSA standards in molecular grade  $H_2O$  and aliquot into two sets of PCR strip tubes; one set is kept at 4C (working set), and the second set is stored at -20C.

15 Add **□200** µl of Reagent B to each well



Wrap plate in foil and put on an orbital shaked at 300 rpm for ( 00:15:00 .

- 17 Read 595 nm absorbance on a suitable plate reader.
- 18 Calculate sample concentration based on the standard curve and adjusted for the applied dilution; remember to subtract the lysis buffer only
- 19 Calculate the amount of sample needed for loading of the required number of wells, taking into account the required mixing with 4X LDS loading buffer and 10X Reducing Agent (RA). 10 μg sample per well in 20 μl (for midi gels) is a good starting point for most applications, but this may require adjusting on a case-by-case basis.

Samples	Protein concentration mg/ml (equivalent to μg/μl)	V for loading of 4 gels with 5 μg per well (prepare for 20 μg in 80 μl, load 20 μl per well)	4X LDS	10X RA	LB_top_up	Final V
G7_25	1.549019608	12.9	20	8	39.1	80
G7_26	1.276348039	15.7	20	8	36.3	80
G7_27	0.580882353	34.4	20	8	17.6	80
G7_28	1.21752451	16.4	20	8	35.6	80
G7_29	1.862132353	10.7	20	8	41.3	80
G7_30	1.144607843	17.5	20	8	34.5	80
G7_31	0.387254902	51.6	20	8	0.4	80
G7_32	1.585784314	12.6	20	8	39.4	80
G7_33	1.441176471	13.9	20	8	38.1	80

Example table of calculations for loading of 4 wells (3+1 extra to take loss into account). LB, lysis buffer.



Once the dilutions have been prepared, return the original samples to § -20 °C for long-term storage. It is also possible to store diluted samples in LDS/RA at § -20 °C, thaw and heat for loading at a later time.

- Heat samples at § 70 °C for © 00:10:00 (if prepared in PCR strips, use thermocycler). It is better to use a lower temperature to limit potential protein degradation while still achieving denaturation. Once completed, leave the samples to cool off at § Room temperature; pulse-spin to pull-down any condensed liquid.
- 21 In the meantime, prepare running buffer according to manufacturers conditions (use 20X MOPS or 20X MES stock buffers choice will depend on the size of your proteins of interest; use deionised H2O for dilution).
  - For the inner gel chamber, prepare a separate aliquot of running buffer and supplement with NuPAGE Antioxidant ( \$\subseteq 500 \mu I \text{ per } \subseteq 200 \mu I \text{ )}.
- 22 Load each sample using gel loading tips; load  $\frac{10}{9}$  protein ladder at the start and/or end of the gel (or as required). Run the gel at 120V until the dye front has reached the plate bottom (usual run time  $\frac{10}{9}$  02:00:00 ).



Note that each midi tank can take 4 gels, however, loading of the gels becomes more difficult when they face backwards and may result in imprecision. It is best to stik to 12-well gels or less to minimise sample diffusion which takes place the longer the samples sit without external voltage during loading.

## Transfer, blocking and primary antibody incubation

- 23 Prepare blocking buffer: [M]3 Mass/Volume Percent BSA to 1X TBS/T (Tween-20 at 0.1 %).
- 24 Prepare the required materials for transfer: trays with EZ gel staining reagent, trays with double-deionised H<sub>2</sub>O for filter paper wetting, tissue paper, scalpel, forceps, transfer packs, iBlot2 machine.
- When the gel has finished running, dissassemble the casette, cut excess gel material off and transfer the gel using the iBlot2 system according to the manufacturer's instructions (Programme P3 for © 00:07:00 is a good starting option, but may require changing based on protein of interest). When assembling the gel transfer stack, wet the filter paper in double-deionised H<sub>2</sub>O and tap dry on some tissue paper before putting on top of the gel; this will reduce water-copped corrosion and the formation of "green" streaks on the membrane).
- Once transferred, cut the membrane if necessary and decide on a marking system to be able to track the identify of the membrane if multiple membranes are being transferred (e.g. corner triangles cut with scalpel; number identifying gel number). Transfer the membrane to a box with blocking solution covering the entire gel. Put to rock on a shaker for 45 minutes to 1 hour.
- 27 Stain the transferred gels in EZ Blue solution for © 01:00:00, then destain using deionised H<sub>2</sub>O. Both steps should be performed on a shaker. The gels can be imaged (epi-illumination or 700/800 nm fluorescence) several days later if necessary leave on bench at room temperature once destaining completed.
- 28 Prepare dilutions of primary antibody in blocking buffer according to the manufacturer's instructions.
- Once the membrane has finished blocking, remove the blocking solution (can be reused), and add the primary antibody dilution so that it covers the entire blot surface. Incubate at 8 4 °C with rotation: typically overnight, but can be left for longer if necessary.

# Secondary antibody incubation and ECL detection

- Decant used primary antibody into original tube (for reuse) and return to the fride (primary antibody should be supplemented with 0.02 % sodium azide for continued fridge storage).
- 31

Wash the membrane 5X for 5 minutes each in 1X TBS/T buffer.

- 32 Dilute the required secondary antibody in blocking buffer. We typically use 1:10,000 dilutions of either one of the following antibodies: anti-rabbit IgG HRP-linked antibody (CST #7074S), anti-mouse IgG HRP-linked antibody (7076S).
- 33 Incubate the membrane in secondary antibody for  $\bigcirc$  01:00:00 at & Room temperature on a shaker. Make sure that the solution covers the entire surface of the membrane.
- 34 Decant secondary antibody and store in fridge for one more reuse if necessary.
- 35

10m

1h

Wash the membrane 5X for 5 minutes each in 1X TBS/T buffer.

- ${\it 36} \quad \hbox{Prepare ECL solution according to the manufacturer's instructions. Use 2\,ml\,ECL\,mix\,per\,10\,cm\,x\,8\,cm\,membrane.}$
- 37 Place the membrane on a flat surface (e.g. glass plate), making sure that excess washing solution has been drained off, and apply the ECL solution. Ensure that the blot is covered evenly.
- 38 Incubate for © 00:05:00 , then proceed with detection using a chemiluminescence imaging system.



Tip: make sure the imaging settings are written down in a readme.txt file, alongside experimental information (samples, transfer condition, antibody dilutions and lot numbers, blocking conditions)

39 Once detected, the membrane can be stored in TBS/T in the fridge for up to two weeks if required for repeated incubations.

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited