

Strep-Tactin AP conjugate Western Blot

Snehadri Sinha

Abstract

Protocol for western blot to detect strep-tagged constructs, using chromogenic detection with Strep-Tactin AP (Alkaline Phosphatase) conjugate. Based on protocol from iba life sciences.

Citation: Snehadri Sinha Strep-Tactin AP conjugate Western Blot. **protocols.io**

dx.doi.org/10.17504/protocols.io.gznbx5e

Published: 17 Jan 2017

Protocol

Blotting

Step 1.

Run your protein samples on SDS-PAGE as usual

Step 2.

Soak 6 pieces of Whatman blotting paper, nitrocellulose membrane (all a bit bigger than your gel) and your SDS-PAGE gel in transfer buffer for a few minutes

Step 3.

Assemble the blot: 3 layers of Whatman paper, nitrocellulose membrane (to the positive pole of the device; in the Bio-Rad blotting machine to the bottom), SDS-PAGE gel, 3 layers of Whatman paper

Step 4.

Roll out all possible air bubbles with a tube or some other device

Step 5.

Blot, use appropriate settings according to your blotting machine.

Chromogenic detection with

Step 6.

Block the membrane with 20 mL TBS-blocking buffer: incubate 1 h at RT with gentle shaking or O/N at 4 C.

Step 7.

Wash 3 times with 20 mL TBS-Tween buffer (each step: 5 min, RT, gentle shaking)

Step 8.

After the last washing step, add 10 mL TBS-Tween buffer to the membrane

Step 9.

(Optional: add 10 uL Biotin Blocking Buffer: 10 min, RT, shaking - blocks endogenously biotinylated proteins which could stain sensitively)

Step 10.

Add 2.5 uL Strep-Tactin AP conjugate (1:4000). Incubate 60 minutes at RT with gentle shaking.

Step 11.

Wash 2 times with TBS-Tween buffer (each step: 1 min, RT, gentle shaking)

Step 12.

Wash 2 times with TBS-buffer (each step: 1 min, RT, gentle shaking)

Step 13.

Transfer membrane in 20 mL reaction buffer with NBT and BCIP. Use manufacturer's instructions for NBT and BCIP amounts; e.g. for Promega NBT+BCIP, first add NBT (33uL for each 5 mL of reaction buffer), mix, add BCIP (16.5 uL for each 5 mL of reaction buffer), mix again, then use within 1 hour)

Step 14.

Proceed with chromogenic reaction under shaking until optimal signal:background ratio is achieved.

Step 15.

Stop reaction by washing several times with distilled H₂O.

Step 16.

Air-dry the membrane and store it in the dark.