



Jul 19, 2019

DNA extraction protocol [↗](#)

PLOS One

Unai Baroja<sup>1</sup>, Inazio Garin<sup>1</sup>, Joxerra Aihartza<sup>1</sup>, Aitor Arrizabalaga-Escudero<sup>1</sup>, Nerea Vallejo<sup>1</sup>, Miren Aldasoro<sup>1</sup>, Urtzi Goiti<sup>1</sup>

<sup>1</sup>Department of Zoology and Animal Cell Biology, University of the Basque Country, UPV/EHU, 48940 Leioa, Basque Country

**1** Works for me [dx.doi.org/10.17504/protocols.io.2ptgdnn](https://doi.org/10.17504/protocols.io.2ptgdnn)

Unai Baroja

## EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0219265>

## THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Baroja U, Garin I, Aihartza J, Arrizabalaga-Escudero A, Vallejo N, Aldasoro M, Goiti U (2019) Pest consumption in a vineyard system by the lesser horseshoe bat (*Rhinolophus hipposideros*). PLoS ONE 14(7): e0219265. doi: [10.1371/journal.pone.0219265](https://doi.org/10.1371/journal.pone.0219265)

## BEFORE STARTING

(N = number of samples)

Set the heatblock at **65 °C**

Set the oven at **37 °C**

**Heat C1 solution** at **65 °C** in the heatblock (aliquoted in a 1.5 ml tube)

Preparte pipette tips (200 ul and 1000 ul)

Label N Dry Bead Tubes with sample numbers/names

- 1 Add the pellets to sterilised weigh trays and weigh the sample (note down number of pellets + total weight)
- 2 Add the pellets to the labelled **Dry Bead Tubes with Bead Solution** and give a quick spin
- 3 Add **60 µl** of preheated **C1 solution** to each **Dry Bead Tube** and vortex thoroughly (15 sec)
- 4 Heat the tubes for 15 minutes at **65 °C** in the heatblock



Meanwhile add **250 µl C2 solution** to empty 1.5 ml tubes and label them (for step 7)



Meanwhile add **200 µl C3 solution** to empty 1.5 ml tubes and label them (for step 10)

5 **Tissuelyse** the tubes 10 minutes at freq 20 (2x20) in ***Precellys Tissue Homogenizer***.

6 **Centrifuge** the tubes 13,000g for 3 minutes

7 Transfer the supernatant (**450 µl**) to the 1.5 ml tubes containing **C2 solution** and vortex briefly

8 **Incubate** for 5 minutes at **4 °C**



Meanwhile label N spin filters and 2N collection tubes

9 **Centrifuge** the tubes 13,000g for 1 minute

10 Transfer the supernatant (up to **600 µl**) to the 1.5 ml tubes containing **C3 solution** and vortex briefly

11 **Incubate** for 5 minutes at **4 °C**




Meanwhile label N low-binding 1.5 ml tubes with sample names/numbers


12 **Centrifuge** the tubes 13,000g for 1 minute

13 Transfer the supernatant (up to **750 µl**) to the 2 ml Collection Tubes

14 **Shake** the **C4 solution**, add  **1100 µl** (550+550) to the supernatant and vortex briefly

15 **Load**  **650 µl** onto a Spin Filter, centrifuge 8,000g for 1 minute and discard the flowthrough to a 50ml falcon tube

15.1 Repeat the **step 15** until loading all the mix

16 Add  **500 µl** of **C5 solution** and centrifuge at 13,000g for 1 minute. Check if the liquid spun through and discard the flowthrough



Avoid splashing C5 solution onto the Spin Filter!

17 **Centrifuge** again at 13,000g for 1 minute and place the Spin Filter in a clean 2 ml Collection Tube

18 Add 50 µl **C6 solution** to the center of the Spin Filter and incubate at  **37 °C** 15 minutes

19 **Centrifuge** 13,000g for 1 minute and discard the Spin Filter

20 Transfer the extract to 1.5 low bind tubes and distribute the extracts in PCR strips



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited