

# Total digestion of marine particles

Paulina Pinedo

## Abstract

Totally digest marine particles collected on polyethersulfone (PES, Pall Supor) filters.

**Citation:** Paulina Pinedo Total digestion of marine particles. **protocols.io**

[dx.doi.org/10.17504/protocols.io.f3wbqpe](https://dx.doi.org/10.17504/protocols.io.f3wbqpe)

**Published:** 11 Oct 2016

## Materials

- ✓ Optima grade Sulfuric acid by Contributed by users
- ✓ Hydrochloric Acid by Contributed by users
- ✓ Hydrogen Peroxide by Contributed by users
- ✓ Hydrofluoric acid by Contributed by users
- ✓ Nitric acid by Contributed by users

## Protocol

### Step 1.

Place filter into PFA vial, ideally with sample side facing downwards or folded inwards.

### Step 2.

Add 1.5 mL of sulfuric acid to all samples and blanks

### Step 3.

Place lid atop loosely and leave filter to soak for 20 min at room temperature

### Step 4.

Lift lid slightly to add 0.5mL of hydrogen peroxide.

### Step 5.

Gently, swirl vial to mix.

### Step 6.

React for 60 min in a hot plate at 110°C with lid placed on loosely

#### **Step 7.**

Add another 0.5 mL peroxide, replace caps loosely and increase hotplate temperature to 200°C

#### **Step 8.**

If undissolved filter pieces or semi-digested viscous material persist, let it cool down and add 100-200µL aliquots of peroxide to digest this material

#### **Step 9.**

Dry vial contents at 235-250°C

Note: If a small droplet of sulfuric acid remains at this point, let vial cool and suspend contents in a small aliquot of Milli-Q water or 8N nitric acid and re-dry at 235°C until the droplet is removed.

#### **Step 10.**

Resuspend dried samples in 2 mL of a freshly prepared mixture of HNO<sub>3</sub>, HCl, and HF acids (4M each) in milli-Q water.

#### **Step 11.**

Heat for 4 hours at 100-110°C.

#### **Step 12.**

Let vials cool to room temperature, uncap and dry at 100-110°C (overnight).

#### **Step 13.**

Add 2mL of freshly prepared 50% HNO<sub>3</sub>/15% H<sub>2</sub>O<sub>2</sub> (v/v)

#### **Step 14.**

Take vials to dryness on a hotplate at 100-110°C

#### **Step 15.**

Resuspend sample in 200-500 µL of 50% HNO<sub>3</sub>/15% H<sub>2</sub>O<sub>2</sub> (v/v). Cap loosely and heat at 110°C.

#### **Step 16.**

After vigorous bubbling cease, uncap vials and dry at 135°C

#### **Step 17.**

Add 100  $\mu\text{L}$  concentrated  $\text{HNO}_3$  and heat at  $110^\circ\text{C}$  to dryness

**Step 18.**

Redissolved in 1mL of 0.1 M  $\text{HNO}_3$

**Step 19.**

Dilute sample to 5% of original concentration with 0.1 M  $\text{HNO}_3$  for analysis

Isotopic Analysis

**Step 20.**

After concentration analysis, follow Conway et al. (2013) to prepare samples for isotopic analysis.