

Small-scale silencing experiment in vegetative *Euplotes crassus* (provisional) Version 3

Rachele Cesaroni

Abstract

Citation: Rachele Cesaroni Small-scale silencing experiment in vegetative *Euplotes crassus* (provisional). **protocols.io**
dx.doi.org/10.17504/protocols.io.hi2b4ge

Published: 01 Apr 2017

Protocol

Step 1.

Grow RNase III deficient *E. coli* strain HT115 in LB with antibiotic selection overnight at 37°C.

Step 2.

Prepare a 1:100 dilution of the bacterial culture, and grow it at 37°C until it reaches an OD₆₀₀ of 0.4.

Step 3.

Add 0.4 mM IPTG and induce RNA transcription from the L4440 plasmid in the bacteria overnight.

Step 4.

Collect the bacteria by centrifugation at max speed for 10 minutes and wash them twice with ddH₂O.

Step 5.

Isolate at least 30 well-starved *Euplotes crassus* cells and feed them twice a day with feeding bacteria.

Step 6.

Leave the *Euplotes crassus* cells at 24°C and monitor daily the change in phenotype.