

# A-Tailing with Taq Polymerase

## New England Biolabs

### Abstract



This protocol can be used to add As to the blunt-ends of DNA fragments that have been amplified using a high-fidelity polymerase (such as Q5® High Fidelity DNA Polymerase).

**Citation:** New England Biolabs A-Tailing with Taq Polymerase. **protocols.io**

dx.doi.org/10.17504/protocols.io.crvv65

**Published:** 03 Feb 2015

## Materials

-  ThermoPol Reaction Buffer Pack - 6.0 ml [B9004S](#) by [New England Biolabs](#)
-  Taq DNA Polymerase with ThermoPol Buffer - 400 units [M0267S](#) by [New England Biolabs](#)

## Protocol

### Step 1.

Clean-up the amplified DNA from the PCR components

#### NOTES

**New England Biolabs** 03 Feb 2015

This can be done by using a PCR-column purification protocol.

**New England Biolabs** 03 Feb 2015

This step is essential because the robust exonuclease activity associated with the high-fidelity enzyme will remove any untemplated nucleotides that are added by Taq DNA Polymerase.

### Step 2.

Set-up the reaction by adding the following components:

#### PROTOCOL

#### . [A-Tailing with Taq Mixture](#)

CONTACT: [New England Biolabs](#)

#### Step 2.1.

PCR-amplified DNA – X


#### Step 2.2.

10X ThermoPol® Buffer – 5ul

#### AMOUNT

5 µl Additional info:

#### REAGENTS

-  ThermoPol Reaction Buffer Pack - 6.0 ml [B9004S](#) by [New England Biolabs](#)

### Step 2.3.

1 mM dATP – 10 µl

### Step 2.4.

Taq DNA Polymerase - **0.2 µl**

 **AMOUNT**

2 µl Additional info:

 **REAGENTS**

 Taq DNA Polymerase with ThermoPol Buffer - 400 units [M0267S](#) by [New England Biolabs](#)

### Step 2.5.

H<sub>2</sub>O to **50 µl**

 **NOTES**

**New England Biolabs** 03 Feb 2015

This volume can be adjusted based on the volume of PCR-amplified DNA that needs to be added

### Step 3.

Incubate the reaction at 72 °C for 20 minutes

 **DURATION**

00:20:00

 **ANNOTATIONS**

**Chi-Yu Lee** 29 Jan 2018

Dear Protocol Author,

I am a phd student at UCL. I want to use Taq polymerase to add A-tail on DNA fragments and then put these fragments into pCRII vector by topo cloning. Do I need to purify the DNA fragments after adding A-tailing to get rid of taq polymerase and other molecules?

Best Regards

Chi-Yu

**New England Biolabs** 29 Jan 2018

Dear Chi-Yu,

Thank you for your question. Yes, the A-tailing reaction should be cleaned up to remove the DNA Polymerase and dNTPs to avoid any unwanted polymerization in the ligation step.

Should you need additional support, please contact NEB tech support directly at 1-800-NEB-LABS or [info@neb.com](mailto:info@neb.com)

Best,

NEB

## Warnings

The DNA cleanup step is essential because the robust exonuclease activity associated with the high-fidelity enzyme will remove any untemplated nucleotides that are added by Taq DNA Polymerase.