



Jul 24, 2019

Make NGM medium V.3

Cancer Research UK / Wellcome Trust Gurdon Institute media kitchen¹¹Cancer Research UK / Wellcome Trust Gurdon Institute

1

Works for me

dx.doi.org/10.17504/protocols.io.5szg6f6



Cristian Riccio



MATERIALS

NAME ▾

CATALOG # ▾

VENDOR ▾

double distilled water (ddH₂O)

Sodium chloride meets analytical specification of Ph.Eur Fisher Chemical

S/3160/65

Fisher Scientific

Agar

AGA02

Formedium

Bacto™ Peptone

211677

Thermo Fisher Scientific

1 ☐ NGM_agar_5L.xls

2 Measure approx 4.8L double distilled H₂O in 5L bell jar with a magnetic flea.

3 Add bacto peptone and NaCl

4 Stir until all solutes are dissolved.

5 Dispense approx 972ml NGM media per 1L bottle OR 486ml NGM into a 500ml bottle.

6 Dispense 17g of agar and a magnetic flea to each 1L bottle OR 8.5g of agar and a magnetic flea to each 500ml bottle.

7 Label, date and autoclave.

8

NaCl		15g
bacto peptone	12.5g	

Double distilled H2O	4.86L	
Agar	per 1L bottle	17g



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited