



May 02,  
2019

Working

## Isolation of endothelial cells from umbilical vessels [↗](#)

PLOS One

Thomas Mudersbach<sup>1</sup>, Daniel Siuda<sup>1</sup>, Karin Kohlstedt<sup>1</sup>, Ingrid Fleming<sup>1</sup>

<sup>1</sup>Institute for Vascular Signalling, Centre for Molecular Medicine, Goethe University, Frankfurt am Main, Germany

[dx.doi.org/10.17504/protocols.io.ybdfsi6](https://doi.org/10.17504/protocols.io.ybdfsi6)

Daniel Siuda

### EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0216218>

### THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Epigenetic control of the angiotensin-converting enzyme in endothelial cells during inflammation PLOS ONE

### Reagents

1

- Hank's Balanced Salt Solution (HBSS) (H6648-500ml, Sigma-Aldrich)
- Enzymatic digestion with Dispase II (29585300, Roche) (2.4 U/ml) in Ham's F12
- Ham's F12 with 10% fetal calf serum (FCS), penicillin (50 U/ml) and streptomycin (50 µg/ml) for centrifugation
- Endothelial cell culture medium: MCDB131 (10372-019, Gibco) containing 10% FCS, L-glutamine (10mM), penicillin (50 U/ml), streptomycin (50 µg/ml), epidermal growth factor (0.1 ng/ml), basic fibroblast growth factor (1 ng/ml) and endothelial growth supplement with Heparin (0.4%)

### Endothelial cell isolation

- 2
1. umbilical cord is cut at one end (about 1 cm) to identify the arteries and the vein
  2. the vein or artery is then cannulated with a button cannula connected to a three-way stop cock. The cannula is fixed with a cable connector to the cord
  3. The cord is rinsed with HBSS until the HBSS remains clear and no blood is washed out any longer
  4. Repeat steps 1-3 on the other end
  5. Now both ends are cannulated and the vessel is filled with Dispase II solution
  6. The cord is incubated for 00:50:00 min at 37 °C
  7. Discard the Dispase II solution and replace with Ham's F12 containing 10 % FCS
  8. Carefully massage the cord to loosen the endothelial cells
  9. Collect the medium in a suitable centrifugation tube and centrifuge for 00:05:00 at 100 x g
  10. Discard the supernatant and resuspend pellet in endothelial cell culture medium and distribute on a fibronectin-coated cell culture dish. Coating with fibronectin should be performed 00:30:00 min to 01:00:00 min before use. Place the cell culture dish containing your cells in a cell incubator at 37 °C in 5% CO<sub>2</sub>.
  11. After 04:00:00 hr cells are washed with cell culture medium.
  12. 48:00:00 hr to 96:00:00 hr after isolation the cells should have formed a confluent cell layer



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited