



Mar 28,
2019

Working

Home-made TOP10 competent heat-shock cells

Forked from [Home-made TOP10 competent heat-shock cells](#)

Magdalena Julkowska¹

¹King Abdullah University of Science and Technology

[dx.doi.org/10.17504/protocols.io.pexdjfn](https://doi.org/10.17504/protocols.io.pexdjfn)



Magdalena Julkowska
King Abdullah University of Science and Technology



ABSTRACT

The standard protocol to propagate TOP10 Competent cells for Heatshock transformation method

EXTERNAL LINK

openwetware.org

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

SAFETY WARNINGS

BEFORE STARTING

Make TSS buffer:

- 5g PEG 8000
- 1.5 mL 1M MgCl₂ (or 0.30g MgCl₂*6H₂O)
- 2.5 mL DMSO
- Add LB to 50 mL
- Filter sterilize (0.22 µm filter)

Note 1: PEG 3350 can be used instead of PEG 8000. According to the original CT Chung paper, PEG 3350 produces better efficiency, and other sized PEGs can be used as well albeit with a slight loss in efficiency.

Note 2: pH from original Chung CT paper calls for pH to be acidic (pH=6.5)... be sure to pH solution before sterilization.

Note 3: If using non-chemically resistant filters (e.g., cellulose nitrate), add DMSO after sterilization. DMSO should be sterile in and of itself, so it may be prudent to add it afterwards if you are unsure about the compatibility of your filters.

- 1 Grow a 5ml overnight culture of cells in LB media. In the morning, dilute this culture back into 25-50ml of fresh LB media in a 200ml conical flask. You should aim to dilute the overnight culture by at least 1/100. Grow the diluted culture to an OD₆₀₀ of 0.2 - 0.5. (You will get a very small pellet if you grow 25ml to OD₆₀₀ 0.2)
- 2 Put eppendorf tubes on ice now so that they are cold when cells are aliquoted into them later. If your culture is X ml, you will need X tubes. At this point you should also make sure that your TSS is being chilled (it should be stored at 4°C but if you have just made it fresh then put it in an ice bath).
- 3 Split the culture into two 50ml falcon tubes and incubate on ice for 10 min.
- 4 Centrifuge for 10 minutes at 3000 rpm and 4°C.

- 5 Remove supernatant. The cell pellets should be sufficiently solid that you can just pour off the supernatant if you are careful. Pipette out any remaining media.
- 6 Resuspend in chilled TSS buffer. The volume of TSS to use is 10% of the culture volume that you spun down. You may need to vortex gently to fully resuspend the culture, keep an eye out for small cell aggregates even after the pellet is completely off the wall.
- 7 Add 100 μ l aliquots to your chilled eppendorfs and store at -80°C . The original paper [1] suggests freezing the cells immediately using a dry ice bath. I (BC) have used liquid nitrogen quite successfully instead of dry ice. Simply placing the cells at -80°C also seems to work well (Jkm) If you run a control every time you clone (i.e. a vector-only ligation), you can as well freeze cells in 200 μ l aliquots. Unused cells can be frozen back once and reused, albeit with some loss of competence.
- 8 It is a good idea to run a positive control on the cells.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited