



## ChroPlate - ProteinA

Alexandra Ehl<sup>1</sup>, David Frommholz<sup>1</sup>, Nadine Stefanczyk<sup>1</sup>

<sup>1</sup>DALEX Biotech

[dx.doi.org/10.17504/protocols.io.uycexsw](https://doi.org/10.17504/protocols.io.uycexsw)



David Frommholz  
DALEX Biotech



### ABSTRACT

#### **Purification Guide for the Isolation of Antibodies with ChroPlate Filtration Plates by DALEX Biotech.**

Easy and quick high throughput antibody purification from various sources and species.

Each well of the ChroPlate has a binding capacity of > 1 mg (tested with human polyclonal Ig, binding varies between species and clones).

The proprietary resin does not shrink or swell in aqueous buffers.

High pressure stability.

pH stability short term 2 - 8, long term 3 - 8.

Excellent thermal stability up to 15 minutes at 80 °C in aqueous buffers at neutral pH.

Can be dried for long term storage (80 °C for > 2 h).

### EXTERNAL LINK

<https://dalex-biotech.com/>

### PROTOCOL STATUS

#### **Working**

Official product protocol by DALEX Biotech.

### GUIDELINES

For optimal binding and purity, the pH of the sample should be 7.5-8.5 and should contain 150-300 mM NaCl. An easy way to achieve this is by adding 1/11 volume 0.5 M Tris, 2 M NaCl (pH 8.0) to your sample. For screening of binding condition for e.g. a monoclonal antibody this parameters might be varied.

Purification works best with an antibody concentration of 1-2 mg/ml in your sample.

It is advisable that all fractions are collected (Sample, flow through, wash, and eluate) in separate plates for analysis, e.g. SDS-PAGE.

### MATERIALS TEXT

Materials provided in the kit:

ChroPlate

Dummy plate

Wash buffer

Elution buffer

Neutralization buffer

Materials not provided in the kit:

Tween-20

0.5 M Tris, 2 M NaCl (pH 8.0)

### SAFETY WARNINGS

The buffers in the kit include sodium azide (CAS No. 26628-22-8) as a preservative.  
For safety information on this chemical(s) check <http://www.dguv.de/ifa/gestis-database>

#### BEFORE STARTING

Make sure your sample is free of particulate matter. You can remove particles by centrifugation or filtration (0.45 µm).

#### How do you want to purify?

- 1 Do you want to purify by centrifugation or by vacuum filtration? Please choose below.

step case

##### Centrifugation

no description provided



#### Equilibration

- 2 Add 500 µl wash buffer to each well, place the ChroPlate on top of a deep-well plate, and centrifuge 5 minutes at 1000 g in a swing-out rotor. For counterbalance of the centrifuge a dummy filter plate is included in the kit.

#### Load and Wash

- 3 Place the ChroPlate on a clean deep-well plate. Add up to 1.5 ml sample to every well. Centrifuge 5 minutes at 1000 g in a swing-out rotor.

##### NOTE

It is advisable to keep the flow-through and wash fractions for later analysis, e.g. SDS-PAGE.

The centrifugation time depends on the sample's volume and viscosity. Volumes larger than 1 ml and viscous samples e.g. serum usually require more than 5 minutes centrifugation time.

For optimal binding and purity, the pH of the sample should be 7.5 - 8.5 and should contain 150 - 300 mM NaCl. An easy way to achieve this is by adding 1/11 volume 0.5 M Tris, 2 M NaCl (pH 8.0) to your sample.

- 4 Empty the deep-well plate or place the ChroPlate on a clean one. Add 500 µl wash buffer to each well and centrifuge 5 minutes at 1000 g in a swing-out rotor.

##### NOTE

In case of unspecific hydrophobic and/or ionic interactions include up to 1 % Tween-20 and/or up to 0.5 M NaCl in the wash buffer.

- 5 Repeat the previous step.

##### NOTE

For increased purity, repeat the washing step a third time.

#### Elution

- 6 Place the ChroPlate on a clean deep-well plate. Add 100 µl elution buffer to each well and centrifuge 5 minutes at 1000 g in a swing-out rotor. Repeat two more times.

- 7 Add one drop of neutralization solution to every well of the deep-well plate and mix gently.

## Equilibration

step case

### Vacuum Filtration

no description provided

- 2 Assemble your vacuum device according to the manufacturer's instructions.  
Add 500 µl wash buffer to each well, place the ChroPlate on a waste receptacle and apply a vacuum of 100 mbar for 2 minutes.

## Load and Wash

- 3 Place the ChroPlate on a clean deep-well plate. Add up to 1.5 ml sample to every well. Apply a vacuum of 100 mbar for 2 - 3 minutes.

#### NOTE

It is advisable to keep the flow-through and wash fractions for later analysis, e.g. SDS-PAGE.

Strength and time of the vacuum depends on the sample's volume and viscosity. Volumes larger than 1 ml and viscous samples e.g. serum usually require more than 2 minutes.

For optimal binding and purity, the pH of the sample should be 7.5 - 8.5 and should contain 150 - 300 mM NaCl. An easy way to achieve this is by adding 1/11 volume 0.5 M Tris, 2 M NaCl (pH 8.0) to your sample.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited