

Protocol C: Preparation of heparin for transfection

David Booth

Abstract

Citation: David Booth Protocol C: Preparation of heparin for transfection. **protocols.io**

dx.doi.org/10.17504/protocols.io.h8vb9w6

Published: 22 Jun 2018

Protocol

Make a 100 mg/ml solution of heparin in sterile water

Step 1.

- Dissolve 5 g of solid heparin sodium (Sigma H3149-25KU) in 45 ml of sterile water

Sterile filter heparin solution through a 0.22 µm filter

Step 2.

Make 1 ml aliquots and store at -20°C

Step 3.