



Jan 06,
2020

mRNA purification with Ambion Dynabeads

Magdalena Julkowska¹

¹King Abdullah University of Science and Technology

1 Works for me

dx.doi.org/10.17504/protocols.io.pbrdim6

Salt Lab KAUST



Magdalena Julkowska
King Abdullah University of Science and Technology



ABSTRACT

The mRNA isolation that is usually used prior to RNAseq library preparation protocol

Adjust the volume of 2 ug of RNA to 36 uL with distilled water

1

Heat to 65°C for 2 minutes to disrupt secondary structures

2

Transfer 72 uL of well resuspended Dynabeads to microcentrifuge tube - place on magnet rack for 30 seconds or until all Dynabeads migrate to the tube well

3

Discard the supernatant, remove tube from the magnet rack, and add 100 uL Binding Buffer to calibrate the beads. Put back on the magnet rack and remove the supernatant

4

Add 36 uL Binding Buffer to Dynabeads. Optimal hybridization required Binding Buffer added in 1:1 ratio relative to sample volume.

5

Add the total RNA to the Dynabeads / Binding Buffer suspension. Mix thoroughly and rotate on a mixer for 5 minutes at room temperature to allow mRNA to anneal to the oligo dT25 on the beads.

6

Place the tube on the magnet rack and wash mRNA-bead complex twice with 200 uL Washing Buffer B. Remove all the supernatant between each washing step with help of the magnet.

7

Elute the samples with 10-20 uL of 10 mM Tris-HCl pH 7.5. Heat to 65°C for 3 minutes and place the tube immediately on the magnet rack.

8

Transfer the eluted mRNA to a new RNase free tube

9



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited