



Oct 18, 2019

Fluorescence measurement

iGEM Dusseldorf¹¹Heinrich-Heine Universität Düsseldorf

1

Works for me

dx.doi.org/10.17504/protocols.io.8f9htr6



iGEM Dusseldorf ⚡

MATERIALS TEXT

LB-medium

antibiotic resistance

Tergitol

24 well plate

96 well plate

- 1 Prepare 5 ml LB medium with antibiotic resistance and 0.5% ml Tergitol
- 2 Inoculate one colony in the prepared LB medium
- 3 Add 0.01-1 mM of the desired fatty acid to the culture medium
- 4 Inoculate the culture for at least 16 hours at 37 °C
- 5 Measure the absorption with a photometer at a wavelength 600 nm
- 6 Dilute the culture to an OD₆₀₀ of 0.05 with LB medium
- 7 Fill a 24 well plate with 800 µl of 3 biological with 3 technically replicats respectively like shown in the table below

8

Sample 1.1	Sample 1.2	Sample 1.3	Sample 5.1	Sample 5.2	Sample 5.3
Sample 2.1	Sample 2.2	Sample 2.3	Sample 6.1	Sample 6.2	Sample 6.3
Sample 3.1	Sample 3.2	Sample 3.3	Sample 7.1	Sample 7.2	Sample 7.3
Sample 4.1	Sample 4.2	Sample 4.3	Sample 8.1	Sample 8.2	Sample 8.3

- 9 Inoculate the 24 well plate at 37°C for at least 16 hours
- 10 Transfer the samples from the 24 well plate on a 96 well plate like shown in the table below

Sample 1.1	Sample 1.2	Sample 1.3	Sample 2.1	Sample 2.2	Sample 2.3	Sample 3.1	Sample 3.2	Sample 3.3	Sample 4.1	Sample 4.2	Sample 4.3
Sample 1.1	Sample 1.2	Sample 1.3	Sample 2.1	Sample 2.2	Sample 2.3	Sample 3.1	Sample 3.2	Sample 3.3	Sample 4.1	Sample 4.2	Sample 4.3
Sample 1.1	Sample 1.2	Sample 1.3	Sample 2.1	Sample 2.2	Sample 2.3	Sample 3.1	Sample 3.2	Sample 3.3	Sample 4.1	Sample 4.2	Sample 4.3
Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank
Sample 5.1	Sample 5.2	Sample 5.3	Sample 6.1	Sample 6.2	Sample 6.3	Sample 7.1	Sample 7.2	Sample 7.3	Sample 8.1	Sample 8.2	Sample 8.3
Sample 5.1	Sample 5.2	Sample 5.3	Sample 6.1	Sample 6.2	Sample 6.3	Sample 7.1	Sample 7.2	Sample 7.3	Sample 8.1	Sample 8.2	Sample 8.3
Sample 5.1	Sample 5.2	Sample 5.3	Sample 6.1	Sample 6.2	Sample 6.3	Sample 7.1	Sample 7.2	Sample 7.3	Sample 8.1	Sample 8.2	Sample 8.3
Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank	Blank

12 Measure the desired fluorescence at the appropriate wavelength and the OD₆₀₀ with a plate reader



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited