

# Oil Red O Staining Drosophila Larval and Prepupal Tissues

Sonia Hall

## Abstract

**Citation:** Sonia Hall Oil Red O Staining Drosophila Larval and Prepupal Tissues. **protocols.io**  
dx.doi.org/10.17504/protocols.io.euhbet6  
**Published:** 08 Apr 2016

## Protocol

### Step 1.

Dissect tissues in 1X PBS and fix in 4% Paraformaldehyde for 20 minutes

### Step 2.

Rinse twice with 1X PBS

### Step 3.

Stain for 20 minutes in Oil Red O Solution

#### ■ ANNOTATIONS

**Lenny Teytelman** 08 Apr 2016

See the "[Oil Red O Solution](#)" for details.

### Step 4.

Rinse 2X in water

### Step 5.

Mount in glycerol and image immediately