



Jan 06,
2020

Enzymatic Assay of Protease Using Azocasein as Substrate

Neilier Junior¹, Rafael de Almeida Barros¹, Camilo Elber Vital¹, Samuel Lessa Barbosa¹, João Vitor Aguiar de Oliveira¹, Gabriele Corrêa da Rocha¹, Cauê Neves Oliveira¹, João Victor Marques Gonçalves Assis¹

¹Universidade Federal de Viçosa

1 Works for me dx.doi.org/10.17504/protocols.io.bayzifx6



Neilier Junior
Universidade Federal de Viçosa

MATERIALS

NAME	CATALOG #	VENDOR
Calcium chloride	1.02378.0500	Merck Millipore
Trichloroacetic acid (TCA)	T6399	Sigma – Aldrich
Sodium hydroxide	S8045	Sigma – Aldrich
Trizma® base	T4661	Sigma Aldrich
Azocasein	A2765	

SAFETY WARNINGS

Wear personal protective equipment: gloves, lab coat and mask.

BEFORE STARTING

Organize your workspace

Make sure all solutions and equipment are available.

Reagent Preparation

- 100 mM Tris-HCl buffer, pH 8.0, 20 mM CaCl₂, at 37 °C.
 - 2.0% (w/v) Azocasein Solution
Heat gently (do not boil) to 50 - 60 °C for 10 min with stirring.
Adjust the pH to 8.0 at 37 °C, if necessary, with either 1.0 M NaOH or 1.0 M HCl.
 - 110 mM Trichloroacetic Acid Reagent (TCA). Dilute with deionized water.
 - 500 mM Sodium Hydroxide (NaOH) Solution. Prepare in deionized water.

Check how many samples will be analyzed to calculate the required volume of each solution to be prepared.

Procedure

2

Pipette (in microliters) the following reagents into 2.0 mL microtubes.

	Blank	Test
Tris-HCl buffer	750 µL	450 µL
Azocasein	750 µL	750 µL
<i>Mix and equilibrate to the at desired temperature. Then add:</i>	*	
Sample (enzyme source)	-	300 µL
<i>Mix and incubate at desired temperature for exactly 30 min.</i>	*	
<i>Remove a 1 mL aliquot from both (test and blank) solutions and place into 2.0 mL microtubes. Then add:</i>		
TCA	1000 µL	1000 µL
<i>Centrifuge at 20,000 g for 10 min. Remove a 1 mL aliquot from supernatant (test and blank) and place into 2.0 mL microtubes. Then add:</i>	*	
NaOH	1000 µL	1000 µL
<i>Mix and transfer the Test and Blank solutions to suitable cuvettes. Measure the A_{440nm} for Test and Blank using a spectrophotometer.</i>	*	

Calculation

$$3 \quad \Delta A_{440nm} = A_{440nm} \text{Test} - A_{440nm} \text{Blank}$$



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited