

# Cassava lea DNA extraction

Devang Mehta

## Abstract

**Citation:** Devang Mehta Cassava lea DNA extraction. **protocols.io**

[dx.doi.org/10.17504/protocols.io.iatcaen](https://dx.doi.org/10.17504/protocols.io.iatcaen)

**Published:** 03 Jun 2017

## Protocol

### Step 1.

Prepare 1.5ml tube by adding a few glass beads.

### Step 2.

To each tube add 200mg of leaf material

### Step 3.

Freeze in liquid N<sub>2</sub> and grind using the dental grinder

### Step 4.

Add 2% beta-mercaptoethanol to the CTAB buffer just before use (200ul/10ml)

CTAB Buffer:

2%	CTAB
2%	PVP-40
100mM	Tris-HCL pH 8.0
25mM	EDTA-Na
2M	NaCl
0.5g/L	spermidine

Make up to required volume with DEPC treated water.

Incubate for 1h at RT and then autoclave

### Step 5.

Add 1ml CTAB buffer to the sample, mix, incubate 15min at 50°C

### Step 6.

Centrifuge 5min full speed, transfer 900ul supernatant to fresh 2ml tube

### Step 7.

Add 900 ul Chloroform:isoamyl alcohol (24:1), mix, centrifuge for 10min full speed

**Step 8.**

Transfer 800ul supernatant to fresh tube and add equal volume of Chloroform:isoamy alcohol (24:1), mix, centrifuge 10 min full speed. Transfer 650ul supernatant to fresh 1.5 ml tube, add 390 ul (0.6 volume) cold EtOH

**Step 9.**

Incubate at -80°C for 30 min or at -20°C overnight.

**Step 10.**

Centrifuge 30 min full speed at 4°C

**Step 11.**

Remove supernatant

**Step 12.**

Wash the pellet in 1 ml 80% ETOH

**Step 13.**

Centrifuge 5min and vacuum dry the pellet.

**Step 14.**

Dissolve the pellet in 100 ul DEPC-treated water. \*continue on ice

**Step 15.**

Measure the DNA concentration at a NanoDrop.

**Step 16.**

Pipette out as much DNA as required.

**Step 17.**

Continue with the following steps for RNA extraction:20. Bring up the volume of total NA to 300 ul

**Step 18.**

Add 1/5th volume of 10M LiCl

**Step 19.**

Incubate at -20 for at least 1 hour

**Step 20.**

Centrifuge full speed (4C) for 30 min

**Step 21.**

Wash with 80% EtOH

**Step 22.**

Centrifuge 5min at 4000 rpm, at 4C, and vacuum dry the pellet.

**Step 23.**

Resuspend in 30-50ul of H2O