

SQK-MAP006 Low Input protocol for library preparation for Nanopore sequencing

Benjamin Istace, Anne Friedrich, Léo d'Agata, Sébastien Faye, Emilie Payen, Odette Beluche, Claudia Caradec, Sabrina Davidas, Corinne Cruaud, Gianni Liti, Arnaud Lemainque, Stefan Engelen, Patrick Wincker, Joseph Schacherer, Jean-Marc Aury

Abstract

Describes the library preparation for Nanopore sequencing from low input DNA according to the SQK-MAP006 protocol

It accompanies the *GigaScience* publication:

Benjamin Istace, et al. (2017) De novo assembly and population genomic survey of natural yeast isolates with the Oxford Nanopore MinION sequencer. *GigaScience*...

Citation: Benjamin Istace, Anne Friedrich, Léo d'Agata, Sébastien Faye, Emilie Payen, Odette Beluche, Claudia Caradec, Sabrina Davidas, Corinne Cruaud, Gianni Liti, Arnaud Lemainque, Stefan Engelen, Patrick Wincker, Joseph Schacherer, Jean-Marc Aury SQK-MAP006 Low Input protocol for library preparation for Nanopore sequencing. **protocols.io**
dx.doi.org/10.17504/protocols.io.gvwbw7e

Published: 06 Jan 2017

Protocol

DNA fragmentation

Step 1.

Fragment DNA (500ng) to a 8Kb size using Covaris gTube

DNA repair

Step 2.

Perform FFPE treatment (NEBNext® FFPE DNA Repair Mix) of fragmented DNA

Clean Up

Step 3.

The DNA repair reaction was cleaned up with Agencourt AMPure XP beads (1x)

End Repair

Step 4.

100 ng DNA fragments are End-repaired using NEBNext® Ultra™ II End Repair Module

dA-tailing

Step 5.

DNA fragments were dA-tailed by using the NEBNext® dA-Tailing Module

Clean Up

Step 6.

The End-Prep reaction was cleaned up with AMPure beads (1x)

Adaptors ligation

Step 7.

Nanopore adaptors were ligated to the dA-tailed DNA fragment by using NEB Blunt/TA Ligase Master Mix

Clean Up

Step 8.

The ligation reaction was cleaned up with MyOne C1-beads (1x)