



Feb 22, 2019

Working

RNA isolation and cDNA synthesis from the marine ichthyosporean *Sphaeroforma arctica* [↗](#)

PLOS Genetics

Sebastián Najle¹¹Institut de Biologia Evolutiva (CSIC-Universitat Pompeu Fabra)[dx.doi.org/10.17504/protocols.io.wqdfds6](https://doi.org/10.17504/protocols.io.wqdfds6)

Multicellgenomelab



Sebastián Najle

Institut de Biologia Evolutiva (CSIC-Universitat Pompeu Fabr...



ABSTRACT

This protocol describes the steps for isolation of total RNA for cDNA synthesis in the marine Ichthyosporean *Sphaeroforma arctica*. In general, Ichthyosporeans are characterized by having a thick cell wall, of unknown composition, that is particularly hard to break, making cell lysis difficult. Here we have introduced a simple modification on the conventional RNA extraction protocol using TRIzol reagent, by including two cycles of freezing-thawing of the cells-TRIzol suspension, in order to facilitate cell lysis. This modification improved significantly the RNA isolation yield. Total RNA is further treated with DNase I enzyme and purified by precipitation to minimize genomic DNA contamination. Highly pure total RNA is then used to synthesize cDNA for downstream applications.

EXTERNAL LINK

<https://doi.org/10.1371/journal.pgen.1007986>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Ocaña-Pallarès E, Najle SR, Scazzocchio C, Ruiz-Trillo I (2019) Reticulate evolution in eukaryotes: Origin and evolution of the nitrate assimilation pathway. PLoS Genet 15(2): e1007986. doi: [10.1371/journal.pgen.1007986](https://doi.org/10.1371/journal.pgen.1007986)

PROTOCOL STATUS

Working

Collect cells by centrifugation

1

- Remove the cells from the culture flasks with aid of a cell scraper (Nunc).
- Transfer cells to a 15 ml centrifuge tube. Pellet cells by centrifugation at 4500 xg for 5 min at **12 °C**.
- Carefully remove supernatant by pipetting. Be cautious not to perturb the pellet, to avoid losing material.

00:08:00

Lyse cells and extract total RNA

2

- Resuspend the cell pellet with **1 ml** of TRIzol (ThermoFisher Scientific) by gently pipetting up and down, until the pellet is completely dissolved. Incubate suspension for 5 min at room temperature.



Make sure that no pellet remaining is visible at the bottom of the tube after resuspension. See manufacturer's protocol for more details on using this reagent.

- Transfer the TRIzol-suspended material to a new 1.5 ml microcentrifuge tube with safe lock (Eppendorf).
- Freeze the suspension by submerging the tube in liquid nitrogen. Once it is completely freeze (nitrogen stops boiling), transfer the tube to a thermal block pre-warmed at 50 °C and thaw the solution for 5 min.
- Repeat freeze-thawing cycle once.



STOPPING POINT: TRIzol suspensions can be kept at **-20 °C** for up to 2 weeks or at **-80 °C** for long-term storage.

🕒 00:25:00

Separate phases and isolate RNA

- 3
- Add **100 µl** of 1-bromo-3-chloropropane (Sigma-Aldrich). Cap the tube and seal it with parafilm to avoid leaks.



1-Br-3-Cl-propane performs better than *chloroform* for phase separation helping to avoid gDNA contamination (interphase is more 'compact' after centrifugation).

- Shake the tube vigorously by hand for 15 seconds and incubate at room temperature for 15 min.
- Centrifuge at 12,000 xg for 15 min at **4 °C**.
- Transfer the aqueous (superior) phase containing RNA (approximately **600 µl**) to a clean, new 1.5 mL tube.



Interphase contains genomic DNA. So, be careful while pipetting to avoid DNA contamination in the RNA sample. Lower phase contain proteins.

🕒 00:30:00

Precipitate RNA

- 4
- Add **500 µl** of isopropanol. Invert 3-4 times to mix and incubate 15 min at room temperature.
 - Centrifuge at 12,000 xg for 10 min at **4 °C**.
 - Discard supernatant and add 1 ml of 75% ethanol (freshly prepared), vortex briefly and centrifuge at 8,500 xg for 5 min at **4 °C**.
 - Discard supernatant completely and dry at room temperature for about 10 min.



CRITICAL!: Avoid overdrying, as it will make RNA very difficult to resuspend.

🕒 00:40:00

Resuspend RNA

- 5
- Add **30 µl** of nuclease-free water, and resuspend by gently pipetting while heating at **56 °C**. Spin down to collect the whole volume at the bottom of the tube.



Water volume can be adjusted as needed. Typically 30 to 50 µL are used.

🕒 00:05:00

Measure RNA concentration

- 6 Measure RNA concentration by spectrophotometry (for example, using nanodrop instrument).

🕒 00:02:00

DNase I treatment

- 7 Transfer up to 10 µg of total RNA to a clean, new tube. Complete volume with nuclease-free water to reach a final volume of 40.5 µl and then add:

- 5 µl 10X DNaseI buffer (ThermoFisher)
- 2 µl RNaseOUT (ThermoFisher)
- 5 µl DNaseI (ThermoFisher)

Incubate at room temperature for 15 min and inactivate the reaction by adding 1 µl EDTA solution and incubating for 10 min at 🌡 65 °C .

🕒 00:20:00

RNA precipitation

- 8 ▪ Add 5 µl of 5M LiCl plus two volumes (110 µl) of 100% ethanol, previously cooled at 🌡 -20 °C . Precipitate overnight at 🌡 -80 °C .

NEXT DAY:

- Centrifuge at maximum speed for 15 min at 🌡 4 °C . Discard supernatant carefully with pipette.
- Wash once with 1 ml 70% ethanol (cold at 🌡 -20 °C).
- Centrifuge at maximum speed for 5 min at 🌡 4 °C , discard supernatant and dry at room temperature for 10-15 min.



CRITICAL!: Be extremely careful not to overdry the sample. Dry RNA is very difficult to resuspend.

- Resuspend pellet (sometimes not visible) with 20 µl of nuclease-free water, by pipetting up and down. Flick the tube, spin down, and incubate overnight at 🌡 -80 °C for complete resuspension.

🕒 48:00:00 (Includes two overnight incubations.)

Measure RNA concentration

- 9 Measure RNA concentration by spectrophotometry (for example, using nanodrop instrument).

🕒 00:02:00

cDNA synthesis

- 10 ▪ Calculate the appropriate volume of DNase-treated total RNA for 0.5 to 3 µg of RNA. Volume should not exceed 8 µl. Assemble the following reaction preferably in DNA LoBind 0.5 µl tubes (Eppendorf).



If DNA LoBind tubes are not available, use normal PCR tubes.

Total RNA ----- x µl
Nuclease-free water ----- 8 - x µl
oligo dT primer ----- 1 µl
10X dNTP mix ----- 1 µl

- On thermocycler, heat the mix for 5 min at 65°C , then cool down immediately on ice.
- While on ice, add 10 μl of cDNA reaction mix of the following composition:

10X Superscript III buffer ----- 2 μl
25 mM MgCl_2 ----- 4 μl
0.1 M DTT ----- 2 μl
RNaseOUT ----- 1 μl
Superscript III ----- 1 μl

- Homogenize the mixture, quickly spin down and immediately incubate for 75 min at 50°C , followed by 5 min incubation at 85°C to inactivate the enzyme.
- Treat samples with 1 μl RNase H (provided with the SuperScript kit), for 1 hour at 37°C (in a water bath).



Reagents used in this step are from SuperScript III First Strand Synthesis Kit (Invitrogen).

02:30:00



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited