

# Pre-stain protocol for GreenGlo™ Safe DNA Dye, 20,000X in Water

Denville Scientific

## Abstract

GreenGlo™, 20,000X in Water, is a non-carcinogenic and non-toxic alternative to Ethidium bromide used for the detection of nucleic acids in agarose gels. It is as sensitive as Ethidium bromide. There is no toxic DMSO as GreenGlo™ is supplied in water.

GreenGlo™ has fluorescence excitation maxima at 295 nm and 490 nm. The fluorescence emission maxima is similar to EtBr when bound to DNA – at 530 nm.

**Citation:** Denville Scientific Pre-stain protocol for GreenGlo™ Safe DNA Dye, 20,000X in Water. **protocols.io**  
dx.doi.org/10.17504/protocols.io.gp6bvre

**Published:** 12 Dec 2016

## Materials

GreenGlo™ Safe DNA Dye, 20,000X in Water, 500µl [CA3600](#) by [Denville Scientific Inc.](#)

## Protocol

### PRE-STAIN

#### Step 1.

Add GreenGlo™ to melted agarose when the agarose has cooled to 50 to 60°C.

#### 📌 NOTES

**Nicole Clouse** 10 Dec 2016

If post-staining, follow the [Post-stain protocol for GreenGlo™ Safe DNA Dye](#).

### PRE-STAIN

#### Step 2.

Use 4-6 µl per 100 ml of agarose or 2-3 µl per 50 ml of agarose

#### 📌 NOTES

**Nicole Clouse** 09 Dec 2016

IMPORTANT : GreenGlo™ is supplied in 20,000X concentration in water

## DETECTION

### Step 3.

Detect bands under UV illumination (yellow or green gelatin- or cellophane filters is recommended for clearer bands) or non-UV LED illuminators such as Blue Light LED illumination.