

# Establishment of axenic sea-ice diatom cultures, modified from Jaeckisch et al. (2011)

Anders Torstensson, Jody Deming

## Abstract

Table 1. Antibiotic concentrations used for establishing axenic sea-ice diatom cultures.

Antibiotic	Final concentration ( $\mu\text{g/ml}$ )
Ampicillin	50
Gentamycin	15
Streptomycin	125
Chloramphenicol	10
Ciprofloxacin	10
Penicillin	100

## References

**Guillard RRL. 1975.** Culture of phytoplankton for feeding marine invertebrates. In: Smith WL, Chanley MH eds. *Culture of marine invertebrate animals*. New York: Plenum, 29-60.

**Jaeckisch N, Yang I, Wohlrab S, Glöckner G, Kroymann J, Vogel H, Cembella A, John U. 2011.** Comparative Genomic and Transcriptomic Characterization of the Toxigenic Marine Dinoflagellate *Alexandrium ostenfeldii*. *PLoS ONE* **6**(12): e28012.

**Citation:** Anders Torstensson, Jody Deming Establishment of axenic sea-ice diatom cultures, modified from Jaeckisch et al. (2011). **protocols.io**

[dx.doi.org/10.17504/protocols.io.qvzdw76](https://dx.doi.org/10.17504/protocols.io.qvzdw76)

**Published:** 11 Jun 2018

## Protocol

### Step 1.

Pre-treat an optically dense diatom culture with 3 min of vortex-mixing ( $< 8^{\circ}\text{C}$ ) and 30 sec of ultrasonication (80% in 5-s pulses) on ice to remove attached bacteria from diatoms and EPS aggregates.

### Step 2.

Centrifuge cells at  $1000 \times g$  for 10 min ( $< 8^{\circ}\text{C}$ ) and wash the pellet twice with sterile F/2 medium (Guillard, 1975) to remove the majority of loose bacterial cells.

### Step 3.

Add six antibiotics to the medium according to Table 1 (see abstract), and incubate for 5 days at  $-1^{\circ}\text{C}$  and under light ( $20\text{--}45 \mu\text{mol photons m}^{-2} \text{s}^{-1}$ ).

### Step 4.

Centrifuge cells at  $1000 \times g$  for 10 min ( $< 8^{\circ}\text{C}$ ) wash pellet in fresh F/2 medium.

**Step 5.**

Pick < 10 cells to start the axenic cultures (pipetting or plate on F/2 agar to pick single colonies).

---

**Step 6.**

After recovery and growth (2 months), check for bacterial cells using DAPI staining and after agar plating on half-strength BD Difco™ Marine Agar 2216 (< 8°C), kept at seawater salinity using artificial seawater (411 mM NaCl, 9.39 mM KCl, 26.1 mM MgCl<sub>2</sub>·6H<sub>2</sub>O, 28.4 mM MgSO<sub>4</sub>·7H<sub>2</sub>O, 5.01 mM TAPSO buffer).

---

**Step 7.**

If bacteria are still present, repeat protocol. Axenic cultures are usually established within the first or second attempt.

---