# Preparation of Tara Sample DNA From Iron-Chloride Precipitates

#### **Matt Sullivan Lab**

# **Abstract**

Citation: Matt Sullivan Lab Preparation of Tara Sample DNA From Iron-Chloride Precipitates. protocols.io

dx.doi.org/10.17504/protocols.io.c24ygv

Published: 04 Jan 2016

# **Guidelines**

References to 'DNA Extraction of Viruses using Wizard Prep Columns' protocol

#### Needed:

- Filters
- Forceps (bleached)
- Aluminum foil squares
- 50cc tube
- Resuspension buffer (0.1M EDTA 0.2M MgCl2 0.2M Ascorbate Buffer)
- Sterile applicator sticks
- 20 mL filter
- 2x 10mL filter
- 20L seawater
- Parafilm
- · Aluminum foil
- Rotator
- · Rubber bands
- 5mL pipet
- 15mL or 50 mL sterile tube
- Centrifuge @ 1000 rpm, 14rpm, 1,000g and 10,000g
- Stock DNase
- DNase buffer
- EDTA
- EGTA
- Amicon Ultra 100kDa centrifugal concentrator
- Wizard Prep Resin
- Wizard Prep Column
- 3mL Luer Lock Syringe
- 5mL Snap-cap tube
- Plunger
- Syringe barrel
- 80% Isopropanol
- 1.5 mL centrifuge tube

- 50-100 µL TE
- Pico Green assay

# **Materials**

Quant-iT dsDNA Pico Green assay kit (Invitrogen) P7589 by Life Technologies

## **Protocol**

# Samples

#### Step 1.

Locate samples and record inventory number and number of filters (or portions) used

## Samples

# Step 2.

Turn the filters precipitate-side out with bleached forceps that have been rinsed with water and use aluminum foil squares as work surface, discarding after each filter

# Samples

#### Step 3.

Remove filter from 50cc tube and refold so that precipitate side comes in contact with resuspension buffer

#### Samples

# Step 4.

Return filter to tube using forceps and sterile applicator sticks

# Samples

## Step 5.

Keep filters dark and refrigerated until ready to suspend

#### **ANNOTATIONS**

# Bonnie Poulos 11 Jan 2016

There should be moisture in the tube - do not let the filters dry out. Usually when the filters are put into the tube initially, they carry some residual seawater with them that serves this purpose. A milliliter of sterile molecular biology grade water can be added if necessary.

# Resuspension of Iron Chloride Precipitates

# Step 6.

Prepare 1x or 2x resuspension buffer.

#### NOTES

VERVE Team 07 Jul 2015

Check that pH of buffer is pH 6.0-6.5

VERVE Team 07 Jul 2015

Prepare about 20% more than needed

VERVE Team 09 Jul 2015

See 0.1M EDTA - 0.2M MgCl2 - 0.2M Ascorbate Buffer

# Resuspension of Iron Chloride Precipitates

#### Step 7.

Prepare 20mL 1x or 10mL 2x per filter from 20L of seawater

## Resuspension of Iron Chloride Precipitates

# Step 8.

Add resuspension buffer to each filter

# Resuspension of Iron Chloride Precipitates

# Step 9.

Parafilm the tubes and wrap all of them in aluminum foil

# Resuspension of Iron Chloride Precipitates

#### Step 10.

Put foil pack of tubes on rotator, in cold room, using rubber bands to secure, and rotate slowly overnight

**O DURATION** 

15:00:00

## Resuspension of Iron Chloride Precipitates

#### **Step 11.**

To recover resuspended viruses, remove the liquid at the bottom using a 5mL pipet and transfer to a fresh 15mL or 50mL sterile, labeled tube

# Resuspension of Iron Chloride Precipitates

# Step 12.

Using bleached and rinsed forceps or sterile applicator sticks, pull edge of filter up and over lip of tube and secure with the lid

# Resuspension of Iron Chloride Precipitates

# **Step 13.**

Centrifuge 1000rpm, 5min, 18°C to recover liquid left on filter

**O DURATION** 

00:05:00

# Resuspension of Iron Chloride Precipitates

#### **Step 14.**

Remove liquid and add more buffer if filter still has a lot of precipitate clinging to it

#### Resuspension of Iron Chloride Precipitates

#### **Step 15.**

Rotate for several more hours

© DURATION

03:00:00

# Resuspension of Iron Chloride Precipitates

## **Step 16.**

Repeat removal of liquid

#### **DNase Treatment**

#### **Step 17.**

Dilute stock DNase 1:100 in 10x DNase buffer (concentration = 400 U/mL)

#### NOTES

# VERVE Team 14 Jul 2015

See DNase Treatment Protocol for details

#### **DNase Treatment**

# **Step 18.**

Add 1/10th volume of diluted DNase to each sample

#### **P** NOTES

VERVE Team 24 Jun 2015

DNase now at 40 U/mL and in 1x reaction buffer

VERVE Team 14 Jul 2015

See **DNase Treatment Protocol** for details

#### **DNase Treatment**

#### **Step 19.**

Parafilm tubes, wrap them in aluminum foil and attach to tube rotator with rubber bands

#### NOTES

VERVE Team 14 Jul 2015

See DNase Treatment Protocol for details

#### **DNase Treatment**

# Step 20.

Incubate by rotating slowly at room temperature for 2 hours

© DURATION

02:00:00

#### NOTES

VERVE Team 14 Jul 2015

See **DNase Treatment Protocol** for details

#### **DNase Treatment**

# Step 21.

Inactivate DNase by adding EDTA and EGTA to 0.1M final concentration each

## **P** NOTES

VERVE Team 14 Jul 2015

See **DNase Treatment Protocol** for details

#### **DNase Treatment**

# Step 22.

Mix by inverting the tube several times

# **P** NOTES

VERVE Team 14 Jul 2015

See **DNase Treatment Protocol** for details

# Concentration and DNA extraction

# Step 23.

Add the DNase treated and inactivated sample to the top reservoir of Amicon Ultra 100kDA centrifugal concentrators

#### NOTES

VERVE Team 24 Jun 2015

Usually this will be in the 15mL size in 50mL centrifuge tubes

VERVE Team 14 Jul 2015

See <u>DNA Extraction of Viruses using Wizard Prep Column</u> for more details

# **ANNOTATIONS**

levgeniia Prekrasna 21 Oct 2017

Are there other options to concentrate viruses instead of using Amicon Ultra centrifugal concentrators? Is it acceptable here to provide PEG precipitation?

## Concentration and DNA extraction

## Step 24.

Centrifuge the concentrators at 1000 g for 5 minute intervals at 18°C until samples are at less than 2mL each

**O DURATION** 

00:05:00

#### NOTES

VERVE Team 14 Jul 2015

See <u>DNA Extraction of Viruses using Wizard Prep Column</u> for more details

#### Concentration and DNA extraction

### Step 25.

Put 1mL of resin on one Wizard Prep column

#### **■** AMOUNT

1 μl Additional info:

# **P** NOTES

VERVE Team 24 Jun 2015

Use no more than 1mL sample per 1mL of Wizard Prep resin (0.5mL is ideal)

VERVE Team 14 Jul 2015

See DNA Extraction of Viruses using Wizard Prep Column for more details

#### **ANNOTATIONS**

Bonnie Poulos 11 Jan 2016

This should state: use 1ml of resin per 0.5ml of sample for application onto one Wizard Prep column.

# Concentration and DNA extraction

## Step 26.

Thoroughly resuspend resin by shaking vigorously

#### **P** NOTES

VERVE Team 14 Jul 2015

See DNA Extraction of Viruses using Wizard Prep Column for more details

# Concentration and DNA extraction

#### Step 27.

Mix 1mL per 0.5-1.0mL of DNA

#### **P** NOTES

VERVE Team 14 Jul 2015

See <u>DNA Extraction of Viruses using Wizard Prep Column</u> for more details

#### ANNOTATIONS

# Bonnie Poulos 11 Jan 2016

Resin works best at a ratio of 1ml resin to 0.5ml sample. If there is more than 1ml sample per 1ml resin, DNA recovery will be reduced.

# Concentration and DNA extraction

# Step 28.

Add each 1mL of resin to a 3mL luer lock syringe attached to a Wizard column

#### **P** NOTES

## VERVE Team 14 Jul 2015

See <u>DNA Extraction of Viruses using Wizard Prep Column</u> for more details

#### Concentration and DNA extraction

# Step 29.

Push through into a 5mL snap-cap tube

#### NOTES

# VERVE Team 14 Jul 2015

See <u>DNA Extraction of Viruses using Wizard Prep Column</u> for more details

#### Concentration and DNA extraction

# **Step 30.**

Save this tube until DNA quantification.

#### NOTES

# VERVE Team 14 Jul 2015

See DNA Extraction of Viruses using Wizard Prep Column for more details

#### Concentration and DNA extraction

# **Step 31.**

Remove syringe from column

## **P** NOTES

## VERVE Team 14 Jul 2015

See <u>DNA Extraction of Viruses using Wizard Prep Column</u> for more details

#### Concentration and DNA extraction

## **Step 32.**

Remove plunger

#### NOTES

#### VERVE Team 14 Jul 2015

See <u>DNA Extraction of Viruses using Wizard Prep Column</u> for more details

# Concentration and DNA extraction

## Step 33.

Reattach syringe barrel to column

#### NOTES

#### VERVE Team 14 Jul 2015

See <u>DNA Extraction of Viruses using Wizard Prep Column</u> for more details

#### Concentration and DNA extraction

# **Step 34.**

Wash columns with 2mL 80% isopropanol pushing through with the syringe barrel into a waste container

#### NOTES

# VERVE Team 14 Jul 2015

See DNA Extraction of Viruses using Wizard Prep Column for more details

#### Concentration and DNA extraction

# Step 35.

Put columns into 1.5mL centrifuge tube and centrifuge at 10,000 g for 2.5 minutes to remove residual alcohol

**O DURATION** 

00:02:30

NOTES

VERVE Team 14 Jul 2015

See <u>DNA Extraction of Viruses using Wizard Prep Column</u> for more details

#### Concentration and DNA extraction

# **Step 36.**

Discard tube

# **₽** NOTES

VERVE Team 14 Jul 2015

See <u>DNA Extraction of Viruses using Wizard Prep Column</u> for more details

#### Concentration and DNA extraction

# **Step 37.**

Place column into a fresh 1.5mL centrifuge tube and pipet on  $50-100\mu L$  of TE ( $0.02\mu m$  filtered and heated to  $80^{\circ}C$ )

#### NOTES

VERVE Team 14 Jul 2015

See <u>DNA Extraction of Viruses using Wizard Prep Column</u> for more details

#### Concentration and DNA extraction

#### **Step 38.**

Vortex gently (1400rpm) and let sit 1 minute

**O DURATION** 

00:01:00

## **P** NOTES

VERVE Team 14 Jul 2015

See <u>DNA Extraction of Viruses using Wizard Prep Column</u> for more details

#### Concentration and DNA extraction

## Step 39.

Centrifuge at 10,000 g for 1 minute

**O** DURATION

00:01:00

## **P** NOTES

VERVE Team 14 Jul 2015

See <u>DNA Extraction of Viruses using Wizard Prep Column</u> for more details

#### Concentration and DNA extraction

# Step 40.

Transfer extracted DNA to Lo-bind DNA 0.5mL tubes

# **P** NOTES

VERVE Team 14 Jul 2015

See <u>DNA Extraction of Viruses using Wizard Prep Column</u> for more details

#### Concentration and DNA extraction

# Step 41.

Repeat elution above (steps 37-40) one more time.

NOTES

# VERVE Team 14 Jul 2015

See <u>DNA Extraction of Viruses using Wizard Prep Column</u> for more details

# **DNA Quantification**

# Step 42.

Use 1-2 $\mu$ L of the first elution and 4-5 $\mu$ L of the second elution for the <u>Pico Green assay</u> to quantify the extracted DNA

Quant-iT dsDNA Pico Green assay kit (Invitrogen)



# **REAGENTS**

Quant-iT dsDNA Pico Green assay kit (Invitrogen) P7589 by Life Technologies

# **DNA Quantification**

# Step 43.

Use the Excel spreadsheet to calculate the ng/µL DNA for each sample

# **DNA Quantification**

# Step 44.

Store in -80°C ultra-low freezer

# **DNA Quantification**

# Step 45.

Label tubes with date, station, depth and concentration and store in -80°C ultra-low freezer