

# Extraction and Lowry-Assay for determination of *Synechocystis* total protein

Anika Wiegard, Sebastian Triesch

## Abstract

A quick Lowry-Assay for the extraction and determination of total protein from *Synechocystis*.

*Photo credit: Miriam Dreesbach, Institute for Synthetic Microbiology, HHU Düsseldorf*

**Citation:** Anika Wiegard, Sebastian Triesch Extraction and Lowry-Assay for determination of *Synechocystis* total protein. **protocols.io**

dx.doi.org/10.17504/protocols.io.taseiee

**Published:** 07 Sep 2018

## Before start

Remember to prepare BSA standards to absolutely quantify your extracted total protein. A BSA range of 5 µg/ml to 100 µg/ml BSA in water or appropriate media showed good results.

## Materials

- Trichloroacetic acid [View](#) by [P212121](#)
- Sodium Hydroxide S320 by [Thermo Fisher Scientific](#)
- 1kg Sodium Carbonate; Na<sub>2</sub>CO<sub>3</sub> (anhydrous) [RC-126](#) by [G-Biosciences](#)
- Copper (II) sulfate pentahydrate [CDB0063.SIZE.2.5Kg](#) by [Bio Basic Inc.](#)
- Folin & Ciocalteu's phenol reagent [F9252](#) by [Sigma Aldrich](#)
- Potassium sodium tartrate tetrahydrate [S2377](#) by [Sigma Aldrich](#)

## Protocol

### Step 1.

**Sample** 1 ml of your *Synechocystis* culture in 2 ml Eppendorf tubes.

### Step 2.

**Add** 110 µl of 100% trichloroacetic acid and **incubate** the mixture on ice for 20 min.

 **DURATION**

00:20:00 : Incubation on ice

**Step 3.**

**Centrifuge** your mixture for 10 min at 15,000 g at 4 °C.

 **DURATION**

00:10:00 : Centrifugation

**Step 4.**

**Discard** the supernatant thoroughly and **place** the tubes **upright** for 10 min. **Tap** the upright tubes carefully until all liquid is removed.

**Step 5.**

**Resuspend** the pellet carefully in 500 µl of a 1 M NaOH solution. **Vortex** and **incubate** the samples over night (approx. 16 hours) at room temperature.

**Step 6.**

**Prepare** a Lowry solution by **mixing** the following reagents *in the given order*:

- 500 µl K-Na-tartrate (2%)
- 500 µl Cu<sub>2</sub>SO<sub>4</sub>\*5 H<sub>2</sub>O (1%)
- 100 ml Na<sub>2</sub>CO<sub>3</sub> (2%)

**Scale** the total volumes **down** to an appropriate amount. You will need 900 µl Lowry-Mix for each sample. Mix the Lowry-solution on the same day you use it and store it in the fridge in the meantime.

**Step 7.**

In a new 1.5 ml Eppendorf tube, **mix** 100 µl sample (in NaOH) and 900 µl Lowry-Mix.

**Step 8.**

**Add** 100 µl 50% Folin & Ciocalteu's phenol reagent (**diluted** in water). Immediately **incubate in the dark** for 45 min.

 **DURATION**

00:45:00 : Incubation in the dark

**Step 9.**

**Spin down** your incubated solution at 14,000 g for 5 min to remove lipids and cell debris.

 **DURATION**

00:05:00 : Centrifugation

### **Step 10.**

Carefully **transfer** 1 ml of your Sample-Lowry-Folin-Mix into a plastic cuvette and measure extinction at 750 nm.

### **Step 11.**

For absolute quantification, **prepare** a standard BSA concentration range from 5 to 100 µg/ml, and **follow** all steps above with your standard sample.



BSA standard set -> go to step #1