

E1202 Blunting Mixture

New England Biolabs

Abstract

Citation: New England Biolabs E1202 Blunting Mixture. **protocols.io**

dx.doi.org/10.17504/protocols.io.crqv5v

Published: 03 Feb 2015

Protocol

Step 1.

Purified DNA (up to 5 µg): **1-19 µl**

Step 2.

10X Blunting Buffer: **2.5 µl**

 **AMOUNT**

3 µl Additional info:

Step 3.

1 mM dNTP Mix: **2.5 µl**

 **AMOUNT**

3 µl Additional info:

Step 4.

Blunt Enzyme Mix: **1.0 µl**

 **AMOUNT**

1 µl Additional info:

Step 5.

Sterile dH2O up to **25 µl**