

Apr 20, 2019

Working

Preparation of human Red Blood Cells for Scanning Electron Microscopy imaging [↗](#)

PLOS One

Marie Anne Balanant¹¹Queensland University of Technology[dx.doi.org/10.17504/protocols.io.yvhfw36](https://doi.org/10.17504/protocols.io.yvhfw36)Marie Anne Balanant
Queensland University of Technology 

ABSTRACT

SEM imaging can be used to quantify RBC shape. RBCs are sampled from fresh blood samples or from packed RBC units during storage. RBCs are resuspended into PBS to make up 500 μ L of solution at 5% haematocrit. RBCs are then fixed by progressively adding a concentrated 2% glutaraldehyde solution in the cell suspension. The cells are incubated for 30 min at RT and in the dark, before being centrifuged and washed in PBS. This fixation protocol was established in order to limit RBC shape changes due to the presence of glutaraldehyde.

After fixation, a RBC suspension at 2.5% haematocrit is adhered to coverslips coated with poly-D-lysine. A second fixation step is realised on the adhered RBCs by incubating the coverslips in osmium tetroxide in cacodylate buffer (1%) for one hour. The coverslips are then dried by incubating them in an ascending series of ethanol, then by incubating them with hexamethyldisilazane (HMDS) for 30 min. The HMDS incubation step is repeated twice before the coverslips are dried in open air. Finally, the coverslips are gold coated and imaged with a Zeiss Sigma FESEM.




EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0215447>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Geekiyana NM, Balanant MA, Sauret E, Saha S, Flower R, Lim CT, Gu Y (2019) A coarse-grained red blood cell membrane model to study stomatocyte-discocyte-echinocyte morphologies. PLoS ONE 14(4): e0215447. doi: [10.1371/journal.pone.0215447](https://doi.org/10.1371/journal.pone.0215447)

MATERIALS

NAME 	CATALOG # 	VENDOR 
Glutaraldehyde EM Grade 25%	G5882-50ML	Sigma Aldrich
Ethanol 100%		
MilliQ Water		
PBS	17-517Q	Lonza
Poly-D-Lysine	A3890401	Thermo Fisher Scientific
Sodium cacodylate trihydrate	View	ProSciTech
Osmium Tetroxide	View	ProSciTech
Hexamethyldisilazane	View	ProSciTech

MATERIALS TEXT

Fresh or stored red blood cells (RBC)

SAFETY WARNINGS

All fixation steps should be conducted under a fume hood. Check SDS before starting work.

BEFORE STARTING

Ethics clearance needs to be obtained before starting work. Check the laboratory has the required safety equipment and PPE for this protocol.

Coverslip preparation

- 1 Prepare glass coverslips by incubating them in **10 Volume Percent** poly-D-lysine solution in PBS for **00:30:00** at RT
- 2 Rinse twice with milliQ water
- 3 Dry and keep at **4 °C** until use

RBC fixation with glutaraldehyde

- 4 Wash **100 µl** of RBC in **900 µl** PBS (**00:02:00** at **1000g**)
- 5 Resuspend **25 µl** of the RBC suspension in **475 µl** of PBS
- 6 Add **500 µl** of **2 Volume Percent** glutaraldehyde solution in PBS to the cell suspension slowly (drop by drop) to reach a final glutaraldehyde concentration of 1%

⚠ SAFETY INFORMATION

Work under a fume hood when using glutaraldehyde

- 7 Incubate for **00:30:00** in the dark
- 8 Wash with **1 ml** of PBS (**00:02:00** at **1000g**)

Adhesion on coverslip

- 9 Dilute **100 µl** of fixed cells suspension into **400 µl** PBS
- 10 Incubate **200 µl** of diluted cell suspension on a poly-D-lysine coated coverslip (**01:00:00**) in corresponding buffer (in 6-wells plates)
- 11 Rinse slide 2 times with PBS (without scratching/damaging cell layer)

SEM sample preparation

- 12 Incubate in **0.1 Molarity (M)** cacodylate buffer (**00:10:00**) . Repeat step 3 times

⚠ SAFETY INFORMATION

Work under a fume hood during this section, check waste disposal for the different chemicals

- 13 Post-fix with [M]1 Mass Percent osmium tetroxide (⌚01:00:00)
- 14 Rinse twice with milliQ water (⌚00:10:00)
- 15 Dehydrate the sample with ethanol
- Ethanol [M]40 Volume Percent (⌚00:10:00)
 - Ethanol [M]50 Volume Percent (⌚00:10:00)
 - Ethanol [M]60 Volume Percent (⌚00:10:00)
 - Ethanol [M]70 Volume Percent (⌚00:10:00)
 - Ethanol [M]80 Volume Percent (⌚00:10:00)
 - Ethanol [M]90 Volume Percent (⌚00:10:00)
 - Ethanol [M]100 Volume Percent (⌚00:15:00)
 - Ethanol [M]100 Volume Percent (⌚00:15:00)
- 16 Incubate in HDMS (⌚00:30:00). Repeat step twice
- 17 Dry under the hood for ⌚00:30:00 minimum (place filter paper under coverslip to prevent attaching to the culture plate)
- 18 Cover with a layer of gold palladium



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited