



## G6PD Activity [↗](#)

PLOS Genetics

Wen Aw<sup>1</sup>

<sup>1</sup>z3314717@unsw.edu.au

[dx.doi.org/10.17504/protocols.io.rwbd7an](https://doi.org/10.17504/protocols.io.rwbd7an)

Wen Aw

### EXTERNAL LINK

<https://doi.org/10.1371/journal.pgen.1007735>

### THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Aw WC, Towarnicki SG, Melvin RG, Youngson NA, Garvin MR, Hu Y, Nielsen S, Thomas T, Pickford R, Bustamante S, Vila-Sanjurjo A, Smyth GK, Ballard JWO (2018) Genotype to phenotype: Diet-by-mitochondrial DNA haplotype interactions drive metabolic flexibility and organismal fitness. PLoS Genet 14(11): e1007735. doi: [10.1371/journal.pgen.1007735](https://doi.org/10.1371/journal.pgen.1007735)

### PROTOCOL STATUS

**Working**

- 1 Homogenize 2 larvae in 40uL of lysis buffer (20mM Tris-HCL, pH8.0, 3mM magnesium chloride, 1mM EDTA, 0.02% B-mercaptoethanol, Protease inhibitor and 0.1% Triton X100)
- 2 Centrifuge at 16,000g for 15min at 4 degree celcius and collected the supernatant
- 3 Protein concentration is determined by Bradford Assay

### Standard Preparation

- 4 Prepare NADPH standard ranges from 0-0.3mM in buffer
- 5 Buffer containing 5mM of both G6P+ 6PGL, 2mM NADP+
- 6 Buffer containing 5mM of 6PGL, 2mM NADP+

### Assay Running

- 7 Added NADP+ into all well (2.5uL of 80mM +37.5uL of Tris Buffer)
- 8 Added 10uL of sample to make up a final volume of 50uL
- 9 Initiate reaction by adding 50uL of substrate (G6P+6PGL (both 5mM final concentration ) or 6PGL (5mM final concentration)



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited