



Sep 24, 2019

Hornwort RNA extraction

Eftychis Frangedakis¹¹University of Cambridge

1

Works for me

dx.doi.org/10.17504/protocols.io.p28dqhw



Eftychis Frangedakis

University of Cambridge, Plant Sciences, OpenPlant



ABSTRACT

The same protocol used for DNA extraction is used for RNA extraction with the addition of an overnight RNA precipitation step with LiCl.

MATERIALS

NAME	CATALOG #	VENDOR
Polyvinylpyrrolidone	PVP40	Sigma – Aldrich
Chloroform	319988	Sigma

MATERIALS TEXT

Extraction buffer:

100mM Tris-HCL pH 8 (2M) 25ml

1.4M NaCl (5M) 140ml

20mM EDTA pH 8 (0.5M) 20ml

2% CTAB 10gr

0.3% β-mercaptoethanol 150μl/50ml

SAFETY WARNINGS

All solutions used during this protocol were prepared using DEPC water!!

- 1 Grind 0.5-2 g of tissue using mortar and pestle in the presence of liquid nitrogen until finely ground. Transfer frozen ground tissue to a 30 ml tube.
- 2 Add 10 ml of 60 °C extraction buffer and 100 mg PVP-40/g tissue. Mix by inversion and incubate in 60 °C in a water bath for 30 min (or incubate on ice for 15 min).
- 3 Remove from heat, and let cool to room temperature for 4 to 6 min.
- 4 Add 12 ml of chloroform:IAA (24:1) and mix by inversion to form an emulsion.
- 5 After mixing thoroughly, spin at 10,000 rpm for 10 min at RT.
- 6 Transfer aqueous phase to a new 30 ml centrifuge tube. Repeat chloroform:IAA extraction to remove cloudiness (PVP) in aqueous phase. Spin at 7,000 rpm for 10 min at RT.

- 7 Add ½ volume of 5M NaCl to the final aqueous solution recovered. Mix well. Add two volumes of cold (-20 °C) 95% ethanol. Mix by inversion. Place in freezer (-20 °C) for 10 min to accentuate precipitation. Spin at 13,000 rpm for 6 min.
- 8 Optional: Resuspend in 2 ml of TE and repeat step 7.
- 9 Pour off supernatant. And wash pellet with cold (0 to 4 °C) 70% v/v ethanol. Dry pellet.
- 10 Dissolve in 400 µl TE and add 100ul 8M LiCl. Incubate o/n at 4 °C.
- 11 Spin 10 min to pellet RNA . Remove supernatant (contains DNA).
- 12 Dissolve pellet with 400ul DEPC dH2O. Add 40ul NaoAc (pH 4.5) and 1ml 95% EtOH and freeze at -80 oC for 10 min. (optional: wash pellet with 95% EtOH).
- 13 Spin for 10 min at 13K and resuspend pellet in 40ul DEPC dH₂O.
- 14 If necessary (pellet not dissolved): incubate at 60 °C for 10 min and then on ice for 30 min. Pellet debris by spinning tube for 10 min at 13K. Transfer supernatant (RNA) to a new tube. Extraction buffer: 100mM Tris-HCL pH 8 (2M) 25ml.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited