



Jun 26,  
2019

## Biolistic Transformation of *Amphidinium* V.2

Isabel Nimmo<sup>1</sup>, Ellen Nisbet<sup>2</sup>, Adrian Barbrook<sup>2</sup>, Chris Howe.<sup>2</sup>

<sup>1</sup>University of Cambridge, <sup>2</sup>Department of Biochemistry, University of Cambridge

Working

dx.doi.org/10.17504/protocols.io.4r2gv8e



Isabel Nimmo  
University of Cambridge



### ABSTRACT

A method to transform the chloroplast of *Amphidinium carterae* using biolistics.

orcid.org/0000-0002-2469-0611

orcid.org/0000-0002-6975-8640

orcid.org/0000-0003-4487-196X

### MATERIALS TEXT

#### pAmpPSBA

gggcgaattg ggcccgcgt cgcgtgctcc cgccgccat ggccgccgcg ggaattcgat 60  
agattaagaa cgattagatt accttcctta ttattgggaa cgttgaacg atgaattaaa 120  
ctcacagtag gcccttcggg cgccttcga acgacgtatt gacctaccta cctaagtaga 180  
gatcgttggt gcgacgtatc aatacccggt taaacaggtg caactgttgt atggttatcc 240  
tgcaacacct cttagatgct agaattgagtc cataaaggct tggaattctc gcttgtagtt 300  
atctctggtt aggtaatagt gccttcctt cttgagcct cctaccgaaa gtcaattcaa 360  
cagcacagat gagtattatt aagagctcaa acaattgtga tcaaagacga tcgaaggata 420  
cggaagcgta gagaatagac aatcgggaga ctatcattag tctggacaa atgtgtcag 480  
acatgaggtt gaccaggtca ttcattaggt ctccggtcat ttgttccat ctctaccca 540  
gtagagaaaa atccaggtca tatcatagga gatggaactg agagatcgag agaacaaga 600  
cagaaaggag gtgaaaggta gttggcatgt tgatctggcg agtcagaggc atcaaacagg 660  
aacgaaggta atagaatgaa tgataattga tgaaacgaca atgaaagtcc tctaagaag 720  
aggcaactaa taaagatgaa acatatctcg ggggtgcgcc taagtcaggc taagaaagca 780  
tcaaacagat acaaccgctg ggcaggtcaa caagtgatta atcctgatct aaggatgatt 840  
acgacctgag gaaagggtgt tacaatacct catacctgta tgactgatcg gtgctttaat 900  
cacacctgac aagttgcaa ttacattgag taggcattct taatagcttg tagggttcac 960  
accgcagtct aaaggtagaa aggtccaaca atcaggtcat gtgcccaga aaccacgctg 1020  
gtaacctgag acctacctc ccttcgggga tcctatcaa tggctgtgt attcagttcc 1080

cgctttacgg agtgctgtg gcggtggagt ggcaggggta ggcaagtagt cgagtgtga 1140  
gctgtaatgg tgtatccgcg ttaccctcc ttgtacctcg gtgcctcggg aagaagacga 1200  
actctcgta gtgaggccag cgtaaaccct ttccatata tattatgaca agccttattc 1260  
gctctaactc ctgggggtct ttcgtcaaa caatcacttc ttctctaac cgtctttata 1320  
tcggttgggt cggctctctc gtctccac ttctctccct tgctaccga gcttatatca 1380  
ctgctttctt ccttgctct gcagtcgata ttgatggat cctgagcca gttgcagggt 1440  
cccttatcta tggtaacaac atcatctctg gtgctgtat tccttctct aacgcaatcg 1500  
gtgtccattt ctatctctc tgggaatccc ttggtctga tgaatggctc tacaacgggt 1560  
gtacttatca gttcgttggt ttccactct tctcgggtg ctgcgggtgg atgggtcgtg 1620  
aatgggaatt ctctatcgt ctcggtatgc gtccatggat ctctagca ttctccgctc 1680  
caattgcgcg agcagctga gtgttcatta ttaccaat cggtcagggt tccttctccg 1740  
atggtatgcc acttggtatc cagggtactt tcaactcat gctgtcttc caagctgagc 1800  
acaagatct catgatctt ttccacattc ttggttagc aggtgtctc ggtggttctc 1860  
tcttccgc aatgcacggg tctctgtat cctctctct cctgcagag actgctggtt 1920  
ctgagtcct taacaacggg tacgtattcg gtcaagagga tgagactat tccatctctg 1980  
cagcacgcg atatttcggg cgtcttatct tccagtatgc tggctcaac aactccggt 2040  
ctcttactt ctcttgca gctggcctg taattggtat ctggtcacc tctcttggtg 2100  
ttgtactat ggcattcaac cttaacgggt tcaactcaa ccagtctatt ctgacgaat 2160  
ctggtcacta cattaactct tgggtgata tctcaaccg tgctgatctt ggtatcgagg 2220  
taatgcacga gcgtaacgct cataactcc ctctgatct tgcatagaat caacccttt 2280  
atccctctt ttgggactct tagaacgact aggcttctg attaacgatg atgtaatt 2340  
attaagaacg atggctagta ataatgcacg aatcactagt gaattcgcg cgcctgcag 2400  
gtcgaccata tgggagagct cccaacgcgt tggatgata gcttgagtat tctatagtgt 2460  
cacctaaata gctggcgta atcatggtca tagctgttc ctgttgaaa ttgtatccg 2520  
ctcacaattc cacacaacat acgagccgga agcataaagt gtaaagcctg ggggtcctaa 2580  
tgagtgagct aactcacatt aattgcgtg cgctcactgc ccgcttcca gtcgggaaac 2640  
ctgctgtgcc agctgcatta atgaatcggc caacgcgagg ggagaggcgg ttgctgatt 2700  
gggcgtctt ccgctctc gctcactgac tcgctgcgt cggctgtcg gctcggcga 2760  
gcggtatcag ctactcaaa ggcgtaata cggttatcca cagaatcagg ggataacga 2820

ggaaagaaca tgtgagcaaa aggccagcaa aaggccagga accgtaaaaa ggccgcgttg 2880  
ctggcgtttt tccataggt cgcggccct gacgagcatc acaaaaatcg acgtcaagt 2940  
cagaggtggc gaaacccgac aggactataa agataccagg cgtttcccc tggaagctcc 3000  
ctcgtgcgt ctctgttcc gaccctgcc cttaccgat acctgtccgc ctttctcct 3060  
tcgggaagcg tggcgcttc tcatagctca cgctgtaggt atctcagttc ggtgtaggtc 3120  
gttcgtcca agctgggctg tgtgcacgaa ccccccgttc agcccgaccg ctgcgctta 3180  
tcggtaact atcgtctga gtccaacccg gtaagacacg acttatcgcc actggcagca 3240  
gccactgga acaggattag cagagcgagg tatgtaggcg gtgctacaga gttctgaag 3300  
tggtggccta actacggcta cactagaaga acagtattg gtatctgcgc tctgtgaag 3360  
ccagttacct tcggaaaaag agttggtagc tcttgatccg gcaaaaaac caccgtggt 3420  
agcgtggtt ttttgttg caagcagcag attacgcgca gaaaaaagg atctcaagaa 3480  
gatccttga tctttctac ggggtctgac gctcagtga acgaaaactc acgttaaggg 3540  
attttggtca tgagattatc aaaaaggatc ttcacctaga tcttttaaa ttaaaaatga 3600  
agtttaaat caatctaaag tatatatgag taaacttgt ctgacagtta ccaatgctta 3660  
atcagtgagg cacctatctc agcgatctgt ctatttcgt catccatagt tgcctgactc 3720  
cccgtcgtg agataactac gatacgggag ggcttaccat ctggccccag tgcgtcaatg 3780  
ataccgcgag acccagctc accggtcca gatttatcag caataaacca gccagccgga 3840  
agggccgagc gcagaagtgg tctgcaact ttatccgcct ccatccagtc tattaattgt 3900  
tgccgggaag ctagagtaag tagttcgcca gtaaatagtt tgcgcaacgt tgttgccatt 3960  
gctacaggca tcgtggtgc acgctcgtc tttggtatgg cttcattcag ctccggttc 4020  
caacgatcaa ggcgagttac atgatcccc atgttgtga aaaaagcgt tagctcctc 4080  
ggctctccga tcgtgtgac aagtaagtt gccgcagtgt tatcactcat ggttatggca 4140  
gcactgcata attctctac tgtcatgcca tccgtaagat gctttctgt gactggtgag 4200  
tactcaacca agtcattctg agaatagtt atgcggcgac cgagttgctc tgcccgcg 4260  
tcaatacggg ataataccgc gccacatagc agaacttaa aagtgtcat cattggaaaa 4320  
cgttctcgg ggcgaaaact ctcaaggatc ttaccgtgt tgagatccag ttcgatgta 4380  
cccactcgt caccgaactg atcttcagca tctttactt tcaccagct tctgggtga 4440  
gcaaaaacag gaaggcaaaa tgccgcaaaa aagggaataa gggcgacacg gaaatgtga 4500  
atactcatac tcttctttt tcaatattat tgaagcattt atcagggtta ttgtctcatg 4560

agcggataca tatttgaatg tatttagaaa aataa caaaa taggggttcc gcgcacattt 4620

ccccgaaaag tgccacctga tgcggtgtga aataccgcac agatgcgtaa ggagaaaata 4680

ccgcatcagg aaattgtaag cgtaaatatt ttgtaaaat tgcggttaa ttttgttaa 4740

atcagctcat ttttaacca ataggccgaa atcggaacaa tccctataa atcaaaagaa 4800

tagaccgaga taggggtgag tgtgttcca gtttgaaca agagtccact attaaagaac 4860

gtggactcca acgtcaaagg gcgaaaaacc gtctatcagg gcgatggccc actacgtgaa 4920

ccatcacctc aatcaagttt ttggggctcg aggtgccgta aagcactaaa tcggaaccct 4980

aaagggagcc cccgatttag agcttgacgg ggaaagccgg cgaacgtggc gagaaaggaa 5040

gggaagaaag cgaaaggagc gggcgctagg gcgctggcaa gtgtagcggc cacgctgcgc 5100

gtaaccacca caccgcgcgc gcttaatgcg ccgctacagg gcgctccat tcgccattca 5160

ggctgcgcaa ctgttggaag gggcgatcgg tgcgggcctc ttcgctatta cgccagctgg 5220

cgaaaggggg atgtgctgca aggcgattaa gttgggtaac gccaggggtt tcccagtcac 5280

gacgttgtaa aacgacggcc agtgaattgt aatacgactc actata 5326

#### To be completed in advance

- 1 Prepare artificial minicircle (pAmpPSBA) at  $1\text{ mg ml}^{-1}$ . pAmpPSBA contains a bacterial origin of replication and an ampicillin resistance marker, so transformation and growth using a suitable *E. coli* strain, followed by a miniprep, is the simplest method to achieve this. pAmpPSBA also contains a linearised minicircle containing the core region and a modified form of PSBA to confer resistance to the herbicide Atrazine.
- 2 Grow wild type *Amphidinium carterae* in *f*/2 medium to early log-phase.
- 3 Prepare *f*/2 plates (1.5% agarose in *f*/2 medium) and leave to dry for 3-4 hours before storing at  $4^{\circ}\text{C}$ . (This drying step is not strictly necessary, but allows the cells to be transferred to the plates much more quickly, which reduces the subsequent time between plating and shooting).
- 4 Autoclave stopping screens, macrocarriers and 1,550 PSI rupture disks (all supplied by BioRad).

#### Cell Preparation

- 5 Count the cells. For each plate to be transformed (including necessary controls), between  $1\text{-}5 \times 10^7$  cells are needed. Harvest an appropriate volume for the number of plates to be shot.

- 6 Centrifuge the cells at low g (~1,500g) and resuspend the pellet in a very small volume (~500µl) of fresh *f/2* medium.
- 7 Spot the cells onto the centre of the *f/2* plate dropwise. Leave to dry. Cells should form a circle of no more than 2cm diameter to maximise the number of cells in the region bombarded by gold particles. Note that *Amphidinium* does not survive long term on solid plates, so aim to shoot as soon after drying as practical.

#### DNA Preparation

- 8 DNA Precipitation is carried out using Seashell Technology's "DNAdeI™ gold carrier" delivery system. 550nm gold particles are used throughout.
- 9 Sonicate the DNAdeI™ gold particles (supplied at 50mg ml<sup>-1</sup>) to dissociate aggregates. 1-2 minutes should normally be sufficient, but the presence of aggregates can be detected by eye.
- 10 Dilute the gold particles in the supplied 'binding buffer'<sup>\*</sup> to a final concentration of 30mg ml<sup>-1</sup>. Use 0.5mg gold per plate for transformation. (e.g. 30ul gold particles and 20ul binding buffer for 3 plates).
- <sup>\*</sup> Supplied by Seashell Technology with the DNAdeI™ gold carrier kit.
- 11 Add artificial minicircle to gold particles at a ratio of 2.5µg DNA per mg gold. Using very concentrated DNA stock keeps volume low and allows for high concentrations of DNA and gold to maximise DNA precipitation onto the gold particles.

#### DNA Preparation

- 12 Vortex the DNA-gold particle mixture briefly.
- 13 Add an equal volume of 'precipitation buffer'<sup>\*</sup> to the DNA-gold particle mix and vortex briefly. Allow to stand for 3 minutes.
- <sup>\*</sup> Supplied by Seashell Technology with the DNAdeI™ gold carrier kit.
- 14 Centrifuge at 10,000g for 10s to pellet the gold particles.
- 15 Remove the supernatant and add 500µl ice cold 100% ethanol. Vortex briefly.
- 16 Centrifuge at 10,000g for 10s to pellet the gold particles.
- 17 Remove the supernatant. Add 10µl ice cold 100% ethanol for each 0.5mg gold in the preparation. (e.g. add 30µl for the 1.5mg gold prepared for shooting 3 plates).
- 18 Briefly sonicate the solution to minimise aggregation and allow for reliable delivery. Again, 1-2 minutes will normally be sufficient, but sonication should continue until even suspension of the gold particles is achieved.

#### Bombardment Preparation

- 19 Wash the macrocarriers in 70% ethanol and allow to dry.
- 20 Wash each rupture disk and macrocarrier in isopropanol and allow to dry.
- 21 Transfer a 10µl aliquot of the prepared gold particles (see above) to the centre of each macrocarrier. Allow to dry.

#### Bombardment

- 22 *Bombardment is carried out in a Biorad 'Biolistic PDS-1000/He' device using Biorad stopping screens, macrocarriers and rupture disks.*
- 23 Shooting is carried out using 1,550 PSI rupture disks. (These show higher rates of transformation than 1,100 or 1,350 PSI disks.)
- 24 Shooting is carried out in vacuum at  $\geq 25$  in. Hg.
- 25 Plates for transformation are placed in the middle slot in the Biolistics device (*i.e. the highest possible position for a plate. Higher slots are occupied by the rupture disk and macrocarrier.*)

#### Post Shooting

- 26 *Cells should be resuspended as soon as possible after shooting to minimise cell death.*
- 27 Resuspend all cells from each plate in 5ml fresh *f*/2 medium. As many cells as possible should be washed off the surface of the plate and collected in the desired vessel for growth.
- 28 Leave to recover overnight.
- 29 Dilute the cells as and apply selection criteria. (*e.g. each plate's cells are diluted to a final volume of 50ml with 2µg ml<sup>-1</sup> Atrazine.*)
- 30 Monitor growth by microscopy every 1-2 days. Untransformed *Amphidinium carterae* are typically dead at 10-12 days post-application of Atrazine.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited