

SN Maintenance Medium for *Synechococcus*

JB Waterbury & JM Willey

Abstract

Reference:

JB Waterbury & JM Willey. Isolation and growth of marine planktonic cyanobacteria. 1988. Methods in Enzymology 167: 100-105.

Citation: JB Waterbury & JM Willey SN Maintenance Medium for *Synechococcus*. protocols.io

dx.doi.org/10.17504/protocols.io.c9hz35

Published: 04 Jan 2016

Guidelines

Reference:

JB Waterbury & JM Willey. Isolation and growth of marine planktonic cyanobacteria. 1988. Methods in Enzymology 167: 100-105.

Formulation:

Stock Solution	Amount	Final Concentration
EDTA, disodium salt, 1g/L	5.57 ml	15 µM
K ₂ HPO ₄ , anhydrous, 6.1g/L	2.57 ml	90 µM
NaNO ₃ , 30%	2.57 ml	9 mM
Na ₂ CO ₃ ·H ₂ O, 4g/L	2.65 ml	100 µM
CTMM ¹	1.00 ml	
Vit. B12, 10mg/ml	0.10 ml	1 µg/L
Seawater	750 ml	75%
DDH ₂ O	250 ml	25%

¹ CTMM contains 1.17 µM Na₂EDTA, 8 nM Zn, 5 nM Co, 90 nM Mn, 3 nM Mo, 10 nM Se, 10 nM Ni and 1.17 µM Fe

Note:

Some *Synechococcus* do not grow well in nitrate containing medium and NH₄⁺ can be substituted for NO₃⁻ as a source of nitrogen (try using 250-500 µM NH₄Cl)

Protocol

Step 1.

Mix seawater and double distilled water

Step 2.

Autoclave 30-35 min.

DURATION

00:35:00

Step 3.

Prepare stock solution following formulation in guidelines

Step 4.

Prepare stock solutions and 0.2 µm filter them

Step 5.

Aseptically add each stock solution to seawater, mixing after each addition

Step 6.

Store in the dark

NOTES

VERVE Team 01 Jul 2015

Vit. B12 degrades in the light.