



Jun 25,
2019

Media and reagents for *Seminavis robusta* cultivation and experiments

Lev Tsylin¹, Aaron Turkewitz²

¹California Institute of Technology, ²University of Chicago

Working

dx.doi.org/10.17504/protocols.io.33xgqpn

Protist Research to Optimize Tools in Genetics (PROT-G)



Lev Tsylin



ABSTRACT

This protocol describes the culturing conditions used to grow *S. robusta* in the lab for all experiments.

GUIDELINES

It matters whether you use NanoPure/MilliQ water or glass distilled water. In glass distilled water-based media, *S. robusta* tends to disadhere from the bottom of the culture flask and float more readily than in MilliQ-based media.

MATERIALS

NAME	CATALOG #	VENDOR
Bioactif artificial sea salt	10304	Tropic Marin
Guillard's (F/2) Marine Water Enrichment Solution	G9903-500ML	Sigma Aldrich
Sodium bicarbonate	S6014	Sigma Aldrich
Penicillin G sodium salt	P3032-25MU	Sigma Aldrich
Gentamycin sulfate	345814-1GM-M	Sigma Aldrich
Ampicillin sodium salt	BP1760-25	Fisher Scientific
Streptomycin sulfate	100556	MP Biologicals
Agar	BP1423-500	Fisher Scientific
175 cm ² flasks with sealing caps	156502	Thermo Scientific
25 cm ² flasks with sealing caps	CC7682-4325	USA Scientific
Stericup	SCGPT05RE	Millipore
Cell Scraper	08-100-241	Fisher Scientific
D-sorbitol	M4125	Sigma Aldrich
Mannitol	S1876	Sigma Aldrich

Preparation of F/2 medium

- 1 Dissolve **34.5 g Tropic-Marin Bioactif (TMB)** in **1 L MilliQ water**

- 2 Add **80 mg sodium bicarbonate**
- 3 Autoclave
- 4 After cooling, add **20 ml 50x Guillard's supplement**
- 5 Filter using Stericup into pre-sterilized bottles
- 6 Store at 4C in the dark

Maintenance antibiotics (200x stock)

- 7 Dissolve **1 g penicillin G sodium salt**, **1 g ampicillin sodium salt**, **0.2 g streptomycin sulfate**, and **0.1 g gentamycin sulfate** in a final volume of **20 ml F/2 medium**.
- 8 Aliquot and store at -20 C.

Agar plates

- 9 Dissolve **17.25 g TMB** in **1 L MilliQ**
- 10 Add **80 mg sodium bicarbonate**
- 11 Add **0 µl** **15 g agar**, dissolve while stirring and heating (may require boiling)
- 12 Autoclave
- 13 After cooling below 60 C, add **20 ml 50x Guillard's supplement**, as well as selective antibiotics, if desired.
- 14 Pour approx. **19 ml** per plate

15 Store in the dark at 4 C

Electroporation buffers

16 [M]1 Molarity (M) sorbitol in MilliQ, [M]0.77 Molarity (M) mannitol in MilliQ, or [M]0.77 Molarity (M) mannitol with [M]8 Volume Percent F/2 in MilliQ.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited