

# Immunohistochemistry - Drosophila Embryo

Sonia Hall

## Abstract

**Citation:** Sonia Hall Immunohistochemistry - Drosophila Embryo. **protocols.io**

dx.doi.org/10.17504/protocols.io.cutwwm

**Published:** 29 Jul 2015

## Guidelines

The following is a list of antibody concentrations for use with the protocol 'Immunohistochemistry - Drosophila Embryo'

Custom antibodies were provided by the published authors.

### DSHB Antibodies

| Antigen              | Host Species | Product ID                      | Concentration | Product Link  |
|----------------------|--------------|---------------------------------|---------------|---|
| DE-Cadherin          | rat          | DCAD2                           | 1:16          | <a href="http://dshb.biology.uiowa.edu/cadherin-DE-">http://dshb.biology.uiowa.edu/cadherin-DE-</a>   |
| Coracle              | mouse        | 50:50 mix of C566.9 and C615.16 | 1:50          | <a href="http://dshb.biology.uiowa.edu/Coracle">http://dshb.biology.uiowa.edu/Coracle</a> and <a href="http://dshb.biology.uiowa.edu/Coracle_2">http://dshb.biology.uiowa.edu/Coracle_2</a> |
| Engrailed            | mouse        | 4D9                             | 1:10          | <a href="http://dshb.biology.uiowa.edu/4D9-anti-engrailed-injected">http://dshb.biology.uiowa.edu/4D9-anti-engrailed-injected</a>   |
| Alpha-spectrin       | mouse        | 3A9                             | 1:16          | <a href="http://dshb.biology.uiowa.edu/3A9-323-or-M10-2_2">http://dshb.biology.uiowa.edu/3A9-323-or-M10-2_2</a>   |
| Mmp                  | mouse        | 3A6B4                           | 1:20          | <a href="http://dshb.biology.uiowa.edu/3A6B4">http://dshb.biology.uiowa.edu/3A6B4</a>   |
| Hindsight            | mouse        | 1G9                             | 1:20          | <a href="http://dshb.biology.uiowa.edu/hindsight-protein">http://dshb.biology.uiowa.edu/hindsight-protein</a>   |
| Atp-alpha            | mouse        | a5                              | 1:50          | <a href="http://dshb.biology.uiowa.edu/a5">http://dshb.biology.uiowa.edu/a5</a>   |
| Nervana              | mouse        | Nrv5F7                          | 1:50          | <a href="http://dshb.biology.uiowa.edu/Nervana-protein">http://dshb.biology.uiowa.edu/Nervana-protein</a>   |
| Discs Large          | mouse        | 4F3                             | 1:10          | <a href="http://dshb.biology.uiowa.edu/4F3-anti-discs-large">http://dshb.biology.uiowa.edu/4F3-anti-discs-large</a>   |
| Fasciclin III        | mouse        | 7G10                            | 1:100         | <a href="http://dshb.biology.uiowa.edu/7G10-anti-Fasciclin-III_2">http://dshb.biology.uiowa.edu/7G10-anti-Fasciclin-III_2</a>   |
| Crumbs (concentrate) | mouse        | Cq4-c                           | 1:50          | <a href="http://dshb.biology.uiowa.edu/Crumbs">http://dshb.biology.uiowa.edu/Crumbs</a>   |
| Armadillo            | mouse        | N2 7A1                          | 1:50          | <a href="http://dshb.biology.uiowa.edu/N2-7A1-ARMADILLO">http://dshb.biology.uiowa.edu/N2-7A1-ARMADILLO</a>   |
| Neuroglian           | mouse        | BP 104                          | 1:25          | <a href="http://dshb.biology.uiowa.edu/neuroglian">http://dshb.biology.uiowa.edu/neuroglian</a>   |
| Flamingo             | mouse        | Flamingo #74                    | 1:20          | <a href="http://dshb.biology.uiowa.edu/Flamingo-74">http://dshb.biology.uiowa.edu/Flamingo-74</a>   |
| N-Cadherin           | rat          | DN-Ex #8                        | 1:20          | <a href="http://dshb.biology.uiowa.edu/cadherin-DN-">http://dshb.biology.uiowa.edu/cadherin-DN-</a>   |
| Rho1                 | mouse        | p1D9                            | 1:40          | <a href="http://dshb.biology.uiowa.edu/p1D9-anti-rho1">http://dshb.biology.uiowa.edu/p1D9-anti-rho1</a>   |
| Alpha-catenin        | rat          | DCAT-1                          | 1:10          | <a href="http://dshb.biology.uiowa.edu/catenin-alpha-">http://dshb.biology.uiowa.edu/catenin-alpha-</a>   |
| Draper               | mouse        | Draper 8A1                      | 1:25          | <a href="http://dshb.biology.uiowa.edu/Draper8A1">http://dshb.biology.uiowa.edu/Draper8A1</a>   |
| Delta                | mouse        | C594.9B                         | 1:50          | <a href="http://dshb.biology.uiowa.edu/Delta-extracellular-domain">http://dshb.biology.uiowa.edu/Delta-extracellular-domain</a>   |
| Elav                 | mouse        | ELAV-9F8A9                      | 1:100         | <a href="http://dshb.biology.uiowa.edu/Elav-9F8A9">http://dshb.biology.uiowa.edu/Elav-9F8A9</a>   |

### Custom Antibodies

| Antigen          | Host Species | Product ID | Concentration | Citation - Pubmed ID |
|------------------|--------------|------------|---------------|----------------------|
| Pio-Pio          | rabbit       | Custom     | 1:100         | 12973360             |
| Neuroglian (1B7) | mouse        | Custom     | 1:10          | 2805067              |
| sqh-1P           | guinea pig   | Custom     | 1:250         | 20920606             |
| sqh-2P           | rat          | Custom     | 1:500         | 20920606             |
| Mcr              | guinea pig   | Custom     | 1:400         | 24496625             |
| Contactin        | guinea pig   | Custom     | 1:2000        | 15459097             |
| Kune             | rabbit       | Custom     | 1:1000        | 20407131             |
| Neurexin         | rabbit       | Custom     | 1:500         | 8978610              |
| Uninflatable     | guinea pig   | Custom     | 1:400         | 19818339             |
| Mtf              | guinea pig   | Custom     | 1:500         | 20935638             |
| Squash           | mouse        | Custom     | 1:400         | 20920606             |

## Protocol

### Step 1.

Time collection of embryos to obtain the correct stage desired for analysis.

Collection is completed using apple juice agar plates with a small amount of yeast paste applied to surface.

Our collections are typically conducted at 25 degrees C, unless otherwise specified in published materials and methods.

### Step 2.

Wash embryos from apple juice agar plate and into an embryo wash basket using a small paintbrush and embryo wash solution.

### Step 3.

Place wash basket into a small tray and pour household bleach on top of the embryos until completely covered. Let sit for 3 minutes to remove chorion.

 DURATION

00:03:00

### Step 4.

Rinse bleach from wash basket and embryos. Once bleach is removed, transfer embryos to container with 50:50 4% paraformaldehyde and heptane.

Shake on orbital rotator at 150rpm for 20 mins

 REAGENTS

✓ Heptane CAS 142-82-5 by Contributed by users

 DURATION

00:20:00

### Step 5.

Allow solution to heptane and PFA to separate. Remove bottom phase (PFA) and dispose of as required. Be sure to not disturb or remove the embryos that are trapped in the interface.

Add Methanol in equal volume to the remaining heptane. Shake vigorously for 20 seconds and allow embryos to settle to the bottom of the container.

Remove heptane from the top layer and most of the methanol. Add methanol and remove two more times.

Transfer embryos in methanol to disposable culture tubes (flint glass) and allow to settle.

Remove methanol without disturbing the embryos at the bottom of the tube.

 REAGENTS

✓ Disposable Culture Tubes 14-958-A by Contributed by users

**Step 6.**

Rinse embryos with 750ul of blocking solution, once.

Allow embryos to settle to the bottom of tube.

Remove blocking solution and rinse twice with 1X PBS.

Remove 1X PBS and add 750ul of blocking solution.

Seal tube with parafilm and rock gently for 30 minutes.

 **DURATION**

00:30:00

**Step 7.**

Remove blocking solution and rinse once with 1X PBS.

Remove 1X PBS and add blocking solution. Add primary antibodies in optimized concentrations. (See Optimized Concentrations for Developmental Studies Hybridoma Bank Antibodies for use with *Drosophila* Embryos).

Wrap tube with parafilm and rock overnight at 4 degrees Celsius.

**Step 8.**

Remove blocking solution with primary antibodies.

Rinse three times with 1X PBS. Remove 1X PBS and add 750ul of blocking solution.

Wrap tube with parafilm and rock at room temperature for 30 minutes.

 **DURATION**

00:30:00

**Step 9.**

Remove blocking solution and rinse once with 1X PBS. Remove 1X PBS and add new blocking solution. Add secondary antibodies at optimized concentrations.

Our studies used concentrations of 1:1000.

Wrap tube with parafilm and foil. Rock gently for 2-3 hours at room temperature.

 **DURATION**

03:00:00

**Step 10.**

Remove blocking solution with secondary antibodies. Rinse three times with 1X PBS. Add 750ul of blocking solution.

Wrap tube with parafilm and foil. Rock for 30 minutes.

 **DURATION**

00:30:00

**Step 11.**

Remove blocking solution and rinse once with 1X PBS.

Carefully remove embryos with a glass pipette that has been rinsed with blocking solution and transfer to a microscope slide.

Using Whatman's paper, carefully wick away excess 1X PBS.

Add 35ul of mounting media and cover with coverslip.

Seal with nail polish if desired.