

Wash Buffer (10mM Tris-HCl pH 7.5, 80% EtOH)

Sarah Hessen-Schmidt

Abstract

Citation: Sarah Hessen-Schmidt Wash Buffer (10mM Tris-HCl pH 7.5, 80% EtOH). [protocols.io](https://doi.org/10.17504/protocols.io.fjmbkk6)

[dx.doi.org/10.17504/protocols.io.fjmbkk6](https://doi.org/10.17504/protocols.io.fjmbkk6)

Published: 10 Aug 2016

Protocol

Step 1.

In an autoclaved bottle add 2.5ml of 1M Tris-HCl pH 7.5 solution

[PROTOCOL](#)


[1M Tris-HCl](#)

CONTACT: [Sarah Hessen-Schmidt](#)

Step 1.1.

Dissolve 30.275g Trizma base in 200ml MilliQ

[REAGENTS](#)

 Trizma Base 93362 by Contributed by users

Step 1.2.

Adjust to pH 7.5 with HCl (35ml)

Step 1.3.

Adjust total volume to 250ml with MilliQ

Step 1.4.

Autoclave

Step 2.

Add 200ml 100% EtOH

[REAGENTS](#)

Ethyl alcohol, Pure 200 proof, for molecular biology [E7023](#) by [Sigma Aldrich](#)

Step 3.

Adjust volume to 250ml by adding 47.5ml MilliQ