



## Protocol for DNA Extraction and Quantitative PCR Detection of *Verticillium dahliae* from Soil [↗](#)

PLOS One

David Wheeler<sup>1</sup>, Jeness Scott<sup>2</sup>, Jeremiah Kam Sung Dung<sup>2</sup>, Dennis Allen Johnson<sup>1</sup>

<sup>1</sup>Washington State University, <sup>2</sup>Oregon State University

[dx.doi.org/10.17504/protocols.io.vwie7ce](https://doi.org/10.17504/protocols.io.vwie7ce)

David Wheeler  
Washington State University

### ABSTRACT

This protocol describes the extraction of DNA from relatively large volumes of soil and subsequent real-time quantitative PCR (qPCR). The qPCR protocol is a modification of the protocol presented by :

Wei, F., Fan, R., Dong, H.-T., Shang, W.-J., Xu, X.-M., Zhu, H.-Q., Yang, J.-R., and Hu, X.-P. 2015. Threshold microsclerotial inoculum for cotton *Verticillium* wilt determined through wet-sieving and real-time quantitative PCR. *Phytopathology* 105:220-229.

### EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0211508>

### PROTOCOL STATUS

**Working**

### MATERIALS

NAME	CATALOG #	VENDOR
Non-Fat Powdered Milk	NB0669.SIZE.250g	Bio Basic Inc.
Ethanol (95 - 100%), molecular grade		
♣ DNA Extraction Kit: MP FastDNA™ SPIN Kit for Soil	116560200	
lysing matrix E	116914050	
OneStep™ PCR Inhibitor Removal Kit	D6030	
SYBR Select Mastermix	4472908	

### MATERIALS TEXT

- sieve (No. 10/2.0 mM)
- sealable plastic bags
- 15 or 50 ml conical polyethylene tubes
- KymWipes
- -80°C Freezer
- -20°C Freezer
- lyophilizer
- FastPrep-24 5G (MP Biomedicals; cat no.: 116005500 )
- Centrifuge
- 2.0 ml microcentrifuge tubes
- Applied Biosystems StepOnePlus Real-time PCR System; cat no.:4376600
- Forward Primer: Vd-F929-947: 5'-CGTTTCCCGTTACTCTTCT-3'
- Reverse Primer: Vd-R1076-1094: 5'-GGATTTCGGCCAGAACT-3'
- Qubit Fluorometer (ThermoFisher Scientific; cat no.:Q33226)

#### Pre-Extraction Sample Preparation:

- 1 The soil cores are combined, air-dried and processed through a sieve (No. 10/2.0 mM). Dried soil samples are stored in a cool dark place until use. When stored in the fridge, keep them in sealed plastic bags to prevent the soil from absorbing moisture.

Place a 5-10 g subsample of soil into a 15 or 50 ml conical polyethylene tube. Cover the tube with a KymWipe to allow moisture and air to

- 2 move through the sample and keep the sample in the tube. Freeze the subsample of soil at -80°C for at least 3 hours.
- 3 The frozen subsample is then freeze dried under a vacuum using a lyophilizer until dry. Drying time varies and ranges from 4 to 24 hours, but could be longer depending on soil type.

#### DNA Extraction

- 4 DNA Extraction Kit: MP FastDNA™ SPIN Kit for Soil (MP Biomedicals; cat no.: 116560200; cost is \$230-290 for 50 reactions or 25 soil samples).
  - 100 ml of 100% ethanol (*not included in the kit*) needs to be added to the concentrated SEW-M solution prior to use.
  - Twice as much lysing matrix E (MP Biomedicals; cat no.: 116914050) is used as compared to other kit components.
- 5 Add 500 mg of freeze-dried soil to a lysing matrix E tube (supplied with the kit).a. Use two tubes per soil sample so that DNA is extracted from a total of 1 gram of soil per sample.
- 6 Pre-grind the samples in a FastPrep-24 5G (MP Biomedicals; cat no.: 116005500 ) for 40 sec at a speed setting of 6.0.
- 7 Add 828 µl of sodium phosphate buffer (supplied with the kit) to each tube.
- 8 Add 150 µl of nonfat powdered milk suspension (133 mg/ml) to each tube (this is not supplied with the kit).
- 9 Add 122 µl of MT buffer (supplied with the kit) to each tube.
- 10 Homogenize in the FastPrep machine for 40 seconds at a speed setting of 6.
- 11 Centrifuge at 14,000 X g for 10 min.
- 12 Transfer supernatant to a clean 2.0 ml microfuge tubes (these tubes are not supplied with the kit).
- 13 Add 250 µl PPS (Protein Precipitation Solution, supplied with the kit) and mix by inverting the tube 10 times.
- 14 Centrifuge at 14,000 X g for 5 minutes to pellet the precipitate.a. When transferring the supernatant to the 1 ml tube containing binding matrix (steps 11-12) avoid transferring any of the precipitate.
- 15 Shake the binding matrix suspension to re-suspend and add 1.0 ml to each 15 ml tube.a. The two 0.5 g subsamples will be pooled in step 12, so only one 15 ml tube for each soil sample is needed.
- 16 Add the supernatants from step 10 to the 15 ml tube containing binding matrix.
- 17 Invert by hand for 2 min to allow the DNA to bind to the binding matrix. Place the tubes in a rack for 5 minutes to allow the silica matrix to settle in the tube.
- 18 Remove and discard 0.75 ml of the supernatant by pipetting. Be very careful not to disturb the silica matrix.
- 19 Gently re-suspend the binding matrix in the remaining amount of supernatant by pipetting.

- 20 Transfer 600 µl of the mixture to a spin filter and centrifuge for 1 min at 14,000 X g.
- 21 Empty the catch tube and add 500 µl more of the re-suspended silica matrix mixture. Repeat steps 15 and 16 until all of the silica matrix mixture has been added to the spin filter.
- 22 Add 500 µl of prepared SEWS-M and gently re-suspend the pellet using the force of the liquid from the pipet tip.a. Ensure that the ethanol has been added to the SEWS-M before use.
- 23 Centrifuge at 14,000 X g for 1 minute. Empty the catch tube.
- 24 Add 500 µl of prepared SEWS-M.a. It is not necessary to re-suspend the pellet during this rinse.
- 25 Centrifuge at 14,000 X g for 1 minute. Empty catch tube.
- 26 Without the addition of any liquid, centrifuge again at 14,000 X g for 2 minutes. This additional spin removes any excess alcohol that might be present. Discard the catch tube and place the spin filter in a clean catch tube (provided in the kit).
- 27 Air-dry the spin filter for 5 minutes at room temperature to allow residual ethanol to evaporate.
- 28 Gently re-suspend the binding matrix in 100 µL DES. (DES is the molecular-grade water that is provided with the kit). Incubate at 55°C for 5 min in a heat block.
- 29 Centrifuge at 14,000 X g for 1 min.

#### Removal of PCR Inhibitors

- 30 Post-DNA Isolation Inhibitor Removal: OneStep™ PCR Inhibitor Removal Kit (Zymo Research; cat no.: D6030; \$102 for 50 reactions or 50 soil samples).
  - *Note: The Zymo columns used for removal of the inhibitors contain a wet matrix. Consequently, it may be necessary to spin briefly the columns for 2-3 sec before use if the matrix is adhered to the cap or if most of it is in the upper portion of the tube.*
- 31 The columns are prepped according to the manufacturer's instructions: 1) snap off the base; 2) remove the green cap; 3) insert into a collection tube; and 4) spin in a microcentrifuge at exactly 8,000 X g for 3 min.a. If the HRC matrix is dry add 400-600 µl water prior to prepping the column.
- 32 Move the prepped column to a 1.5 ml centrifuge tube.
- 33 Add the DNA suspension from step 25 to the prepped sample column.
- 34 Centrifuge at 8,000 X g for 1 min.
- 35 Store the DNA in the freezer (-20°C).

#### Quantitative PCR (adapted from Phytopathology 105:220-229)

- 36 **PCR cocktail:** All samples are analyzed in triplicate (a total of three 20 µl reactions for each soil sample). Data are averaged across the


three reactions.

Reagent	[Initial]	Volume/ 20 µl reaction	[Final]
ABI SYBR Select MasterMix	2X	10 µl	1X
Forward Primer: Vd-F929-947	10 µM	2 µl	1 µM
Reverse Primer: Vd-R1076-1094	10 µM	2 µl	1 µM
Molecular-grade water	NA	4 µl	NA
DNA sample	NA	2 µl	NA

ABI SYBR Select Mastermix (Thermo Fisher Scientific; cat no.:4472908)  
 Forward Primer: Vd-F929-947: 5'-CGTTTCCCGTTACTCTTCT-3'  
 Reverse Primer: Vd-R1076-1094: 5'-GGATTTCGGCCAGAAACT-3'

- 37 Standard Curve:** A standard curve should be included with every qPCR run.
1. Extract *V. dahliae* genomic DNA from mycelia or conidia grown in pure culture.
  2. Quantify DNA with the Qubit Fluorometer (ThermoFisher Scientific; cat no.:Q33226) according to the manufacturer's instructions.
    - *Another fluorometric quantification method can be used, but UV-absorbance-based quantification should not be used.*
  1. Prepare a ten-fold serial dilution of DNA to use as the template for the standard curve.. The DNA concentration should range from 0.5 ng DNA/µl to 0.5 fg DNA/µl
  2. A no template (water) control should also be included with every qPCR run.

- 38 Thermal Cycler Conditions** (Applied Biosystems StepOnePlus Real-time PCR System; cat no.:4376600)
1. 10 min @ 95°C
  2. 40 cycles of 30 s @ 95°C, 30 s @ 60°C, and 30 s @ 72°C (measure data).
  3. Melt curve analysis from 60°C, gradually increasing 0.3°C/s to 95°C.

 This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited