

Acrolein + PFA perfusion for immunocytochemistry

Cleyde Helena

Abstract

Appropriate tissue fixation is essential for a good quality immunocytochemistry (ICC).

There are several fixation methods, but "whole body" or "target" perfusion of the animal is one of the most efficient methods.

Most perfusing protocols include only paraformaldehyde, but acrolein penetrates tissue more rapidly, improving final fixation of the tissue.

Citation: Cleyde Helena Acrolein + PFA perfusion for immunocytochemistry. **protocols.io**
dx.doi.org/10.17504/protocols.io.tb9eir6

Published: 07 Sep 2018

Protocol

Step 1.

Prepare animal as for Paraformaldehyde fixation; e.g. animal anesthetized with Nembutal, open chest cavity; inject 0.1 cc heparin into the heart apex before first rinse. We clamp the descending aorta to eliminate whole body perfusion. Cut atrium.

Step 2.

Saline nitrite rinse until clear effluent; about 150-200 ml/225-250 g rat.

 [AMOUNT](#)

1000 ml : 0.9% saline

 [AMOUNT](#)

20 g : Sodium Nitrite

Step 3.

Switch to Acrolein/Paraformaldehyde; continue until animal is hard, as judged by the front feet. Requires approximately 150-200 ml/225-250 g rat.

 [AMOUNT](#)

500 ml : 8% PFA

 [AMOUNT](#)

100 ml : 0.5 M Monobasic Potassium Phosphate

 [AMOUNT](#)

100 ml : 0.5M Dibasic Potassium Phosphate

 [AMOUNT](#)

275 ml : distilled water

 [AMOUNT](#)

25 ml : acrolein

SAFETY INFORMATION

add acrolein just before ready to perfuse; if using a sealed bottle or flask, slowly release pressure after doing so and before perfusing

Step 4.

Second saline nitrite rinse; ≈150 ml. No post-fixation necessary.

Warnings

Acrolein is very toxic, similar to tear gas, so wear mask and don't get too close!