

# Removal of Melanin

Jason Stajich

## Abstract

Principle is that CTAB is charging the anionic nucleotides whereby neutral polysaccharides/ melanins are remaining in supernatant. This method also uses urea with the idea that the presence of urea helps to solubilize hydrophobic compounds that would otherwise potentially interact with the hydrophobic core of the CTAB micelles.

**Citation:** Jason Stajich Removal of Melanin. **protocols.io**

dx.doi.org/10.17504/protocols.io.chdt25

**Published:** 02 Oct 2014

## Protocol

### Step 1.

Add water is added to 100- 200  $\mu$ L DNA/RNA solution until a volume of 400  $\mu$ L is reached.

### Step 2.

Add 130 $\mu$ L 5M NaCl.

 **AMOUNT**

130  $\mu$ L Additional info:

### Step 3.

Add 1.6 mL of CTAB-Urea solution

 **PROTOCOL**

. [CTAB-Urea buffer](#)

CONTACT: [Jason Stajich](#)

### Step 3.1.

50 mM Tris-HCl, pH 7.0

### Step 3.2.

1% CTAB

### Step 3.3.

4M Urea

### Step 3.4.

1 mM EDTA

### Step 4.

Mix samples (by hand).

### Step 5.

Incubate overnight at 4°C.

 **DURATION**

15:00:00

**Step 6.**

Centrifuge for 15 minutes at max speed at 4°C.

 **DURATION**

00:15:00

**Step 7.**

Remove the solution. Be very careful in this step!

■ **ANNOTATIONS**

**Jason Stajich** 02 Oct 2014

This is a critical step! Not all DNA will form a pellet at the bottom, but it will also adhere to the side of the tube and some will float on top of the solution. Be very careful in removing the solution!

During this removal it is probably needed to perform additional centrifugation steps to reduce risk of losing DNA.

**Step 8.**

Resuspended in 400 µL 7M guanidine hydrochloride

 **AMOUNT**

400 µl Additional info:

**Step 9.**

Add 2 Vol of EtOH (100%)

**Step 10.**

Incubate on ice for 1 hour

 **DURATION**

01:00:00

**Step 11.**

Centrifuge for 15 min at 4 °C

 **DURATION**

00:15:00

**Step 12.**

Wash with 70% EtOH

**Step 13.**

Wash second time with 70% EtOH

**Step 14.**

Centrifuge 10 min at max speed at RT.

 **DURATION**

00:10:00

**Step 15.**

Remove the supernatant.

**Step 16.**

Dry pellet

**Step 17.**

Resuspend in TE