

0.5M EDTA

Ms Alex Aitken

Abstract

This shows how to make 0.5M EDTA

Citation: Ms Alex Aitken 0.5M EDTA. **protocols.io**

dx.doi.org/10.17504/protocols.io.c5wy7d

Published: 25 Jan 2016

Guidelines

Use vigorous stirring, moderate heat (if desired) and time

Protocol

Step 1.

Mix 18.6 g EDTA in 100ml Milli-Q water

 [AMOUNT](#)

19 g Additional info:

Step 2.

Bring pH to 8.0 using NaOH

 [NOTES](#)

VERVE Team 22 Jun 2015

EDTA will not be soluble until pH reaches 8.0-this will take time (hours)