

NGFR-mNG2 sorting protocol

Benjamin Emert

Abstract

Protocol for preparing NGFR-mNG2 tagged WM989 A6-G3 melanoma cells for sorting.

Citation: Benjamin Emert NGFR-mNG2 sorting protocol. **protocols.io**

[dx.doi.org/10.17504/protocols.io.rh6d39e](https://doi.org/10.17504/protocols.io.rh6d39e)

Published: 05 Jul 2018

Before start

Warm trypsin, media and DPBS.

Materials

7-AAD Viability Staining Solution 420403 by
BioLegend

12x75 mm high clarity polypropylene test
tubes 352063 by Corning

Falcon® 5 mL Round Bottom Polystyrene
Test Tube, with Cell Strainer Snap Cap 352235
by Corning

DPBS no calcium, no magnesium 14190136
by Invitrogen - Thermo Fisher

0.05% Trypsin-EDTA, phenol red 25300054 by
Invitrogen - Thermo Fisher

Propidium iodide staining solution 556463 by
BD Biosciences

DAPI D3571 by Invitrogen - Thermo Fisher

Ultrapure 0.5M EDTA pH 8.0 15575020 by

Invitrogen - Thermo Fisher

✓ TU2% view by Contributed by users

15mL polypropylene centrifuge tubes
352096 by Corning

Protocol

Trypsinize and pellet cells

Step 1.

Wash 5 confluent 10cm plates with 1xDPBS

📌 NOTES

Benjamin Emert 05 Jul 2018

For the NGFR-mNG2 WM989 A6-G3 cell lines you can expect 4-5 million cells on a confluent 10cm dish. Scale up the number of starting plates accordingly.

Trypsinize and pellet cells

Step 2.

Aspirate 1xDPBS and add 2mL 0.05% Trypsin-EDTA. Incubate at 37°C for 3-5 minutes.

📌 AMOUNT

2 ml Additional info:

0.05% Trypsin EDTA

🌡️ TEMPERATURE

37 °C Additional info:

incubate

Trypsinize and pellet cells

Step 3.

Neutralize trypsin with 4mL TU2% media. Transfer dissociated cells to 15mL falcon tube.

📌 AMOUNT

4 ml Additional info: TU2%

media

📌 NOTES

Benjamin Emert 05 Jul 2018

Can use 1 15mL falcon tube for 2.5 10cm plates. Or scale down the volume of trypsin and media to fit more cells into fewer tubes.

Trypsinize and pellet cells

Step 4.

Pellet cells by centrifugation at 1,000 rpm (200rcf) for 3 minutes.

NOTES

Benjamin Emert 05 Jul 2018

During centrifugation, prepare sterile 1xDPBS with 2mM EDTA. Prepare excess.

Benjamin Emert 05 Jul 2018

Alternatively, can centrifuge at 300 rcf

Wash cell pellets

Step 5.

Resuspend cell pellet in 5mL 1xDPBS with 2mM EDTA.

AMOUNT

5 ml Additional info:

1xDPBS with 2mM EDTA

NOTES

Benjamin Emert 05 Jul 2018

Can combine pellets at this point

Wash cell pellets

Step 6.

Pellet cells by centrifugation at 1,000 rpm (200rcf) for 3 minutes.

NOTES

Benjamin Emert 05 Jul 2018

During centrifugation prepare 1xDPBS with 2mM EDTA and live dead staining solution. I use one of the following:

4uL/mL 7-AAD (final concentration 200ng/mL)

2uL/mL Propidium iodide (final concentration 100ng/mL)

1uL/mL DAPI (final concentration 50ng/mL).

Prepare excess for diluting cells if necessary.

Benjamin Emert 05 Jul 2018

During centrifugation prepare sterile 5mL polypropylene FACS tubes with 1mL TU2% for collecting sorted cells.

Stain dead cells

Step 7.

Resuspend cell pellets in 1xDPBS with 2mM EDTA and dead cell stain. I will typically use 6mL for 5-10cm plate's worth of cells.

Filter cells

Step 8.

Remove cell clumps by passing cell suspension through cell strainer cap into sterile 5mL polypropylene

FACS tubes. Place tubes on ice.

Proceed with FACS

Step 9.

Bring with you the following:

- 1xDPBS with 2mM EDTA and dead cell stain for diluting cells
 - 5mL polypropylene FACS tubes with 1mL TU2% for collecting cells
 - Excess FACS tubes and media
-