



Nov 12,
2019

Extraction of gDNA from *Synechocystis* 6803

Dennis Dienst¹

¹Department of Chemistry - Ångström, Uppsala

1 Works for me dx.doi.org/10.17504/protocols.io.k99cz96

CyanoWorld



Dennis DD. Dienst
Department of Chemistry - Ångström, Uppsala



THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

null

MATERIALS

NAME	CATALOG #	VENDOR
RNase A 10mg/ml, DNase and Protease-free	EN0531	Thermo Scientific
Phenol/Chloroform/Isoamyl alcohol (25:24:1)	A156	Carl Roth
2-propanol	6752	Carl Roth

Preparation of buffer solutions

- TE buffer:**
10 mM Tris-HCl
1 mM EDTA ; pH 8.0

TES buffer:
50 mM Tris-HCl
1 mM EDTA, pH 8.0
25% (w/v) Saccharose

Cell harvest & disruption

- Centrifuge 50 mL liquid culture (OD ~ 0.4-1.0) for 7 min at 4000 g and 4 °C
⌚ 00:07:00 4000 g
- Resuspend pellet in 10 mL TE buffer
 - repeat centrifugation (step 1)

⌚ go to step #2
- Repeat wash (step2)

⌚ go to step #3

- 5
- Resuspend pellet in 1 mL TES buffer
 - freeze in liquid nitrogen



wear goggles when working with liquid nitrogen
https://en.wikipedia.org/wiki/Liquid_nitrogen

- 6
- add 5 mg × mL⁻¹ (spatula tip) lysozyme, 3.2 units proteinase K and 100 µL 20% SDS
 - incubate for 1h at 37 °C

🕒 01:00:00

🌡 37 °C

- 7
- incubate for 1h at 60°C under gentle agitation
 - alternatively: 16h (o/n) at 37 °C

🕒 01:00:00

🌡 60 °C

Phenol-chloroform extraction

- 8
- add 1 Vol.
 - add 1 Vol. Phenol/Chloroform/Isoamyl alcohol (25:24:1) and mix
 - centrifuge for 10 min at 12000 g and 4°C

🕒 00:10:00 12000 g

🌡 4 °C

- 9
- transfer upper (aqueous) phase to fresh 1.5 mL tube
 - add 1 Vol. Phenol/Chloroform/Isoamyl alcohol (25:24:1) and mix
 - centrifuge for 10 min at 12000 g and 4°C

🕒 00:10:00 12000 g

🌡 4 °C

- 10
- Repeat step 8
- 🔄 go to step #9

DNA precipitation

- 11
- transfer upper (aqueous) phase to fresh 1.5 mL tube
 - add 0.7 Vol. 2-propanol (isopropanol)
 - incubate for 5 min at room temperature

🕒 00:05:00 room temperature

- 12
- centrifuge for 30 min at 12000 g and 4 °C

🕒 00:30:00 12000 g

🌡 4 °C

- 13
- carefully remove 2-propanol
 - add 300 μ L 70% EtOH
 - air-dry pellet for 1h at room temperature
- (i.e. leave tube open on your bench and enjoy a coffee or two)



DNA pellet might be hardly visible

 **01:00:00** room temperature

- 14
- Resuspend DNA in 30-100 μ L TE buffer

optional: RNA removal

- 15
- add 2 μ L RNase A
 - incubate 30 min at 37 °C
 - incubate 10 min at 72 °C

- 16
- Store DNA at -20 °C



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited