

Dephosphorylation using CIP in Restriction Enzyme Reaction (M0290)

New England Biolabs

Abstract

Protocol for Dephosphorylation of 5'-ends of DNA using CIP in Restriction Enzyme Reaction. Uses the Calf Intestinal Alkaline Phosphatase (CIP - M0290)

Citation: New England Biolabs Dephosphorylation using CIP in Restriction Enzyme Reaction (M0290). [protocols.io](https://doi.org/10.17504/protocols.io.cdhs35)
[dx.doi.org/10.17504/protocols.io.cdhs35](https://doi.org/10.17504/protocols.io.cdhs35)

Published: 03 Oct 2014

Materials

🐛 Alkaline Phosphatase, Calf Intest (CIP) - 1,000 units [M0290S](#) by [New England Biolabs](#)

Protocol

Step 1.

Digest 1–5 µg of plasmid DNA in a **20 µl** reaction as follows:

✓ PROTOCOL

. [Mixture for M0290 CIP](#)

CONTACT: [New England Biolabs](#)

■ ANNOTATIONS

New England Biolabs 25 Jan 2015

Scale larger reaction volumes proportionally.

Step 1.1.

DNA ≥ **1 µl**

■ ANNOTATIONS

Peter Rüthemann 14 Oct 2015

Cas9

Step 1.2.

Restriction Enzyme Buffer (10X) **2 µl**

Step 1.3.

Restriction Endonuclease **1 µl**

Step 1.4.

H₂O, purified to **20 µl**

Step 2.

Incubate at 37°C for 60 minutes or follow manufacturer's recommendations.

 DURATION

01:00:00

Step 3.

Add 1 unit of CIP for every 1 pmol of DNA ends (about 1 µg of a 3 kb plasmid).

■ ANNOTATIONS

Mohd Fahrurrazi Tompong 25 Apr 2018

How much exactly 1 unit of CIP in uliter

Step 4.

Incubate at 37°C for 30–60 minutes.

 DURATION

01:00:00

Step 5.

Purify DNA by gel purification, spin-column or phenol extraction.

Step 6.

Proceed with ligation.