



Cytochrome c oxidase assay [↗](#)

PLOS Genetics

Wen Aw¹

¹z3314717@unsw.edu.au

[dx.doi.org/10.17504/protocols.io.shieb4e](https://doi.org/10.17504/protocols.io.shieb4e)

Wen Aw

EXTERNAL LINK

<https://doi.org/10.1371/journal.pgen.1007735>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Aw WC, Towarnicki SG, Melvin RG, Youngson NA, Garvin MR, Hu Y, Nielsen S, Thomas T, Pickford R, Bustamante S, Vila-Sanjurjo A, Smyth GK, Ballard JW (2018) Genotype to phenotype: Diet-by-mitochondrial DNA haplotype interactions drive metabolic flexibility and organismal fitness. PLoS Genet 14(11): e1007735. doi: [10.1371/journal.pgen.1007735](https://doi.org/10.1371/journal.pgen.1007735)

PROTOCOL STATUS

Working

- 1 Third instar larvae were ground for 10 s in 100 µL of ice-cold homogenization buffer (50 mM phosphate buffer containing 0.05% Tween-80) using a Kontes pellet pestle motor.
- 2 Homogenates were then diluted with an additional 500 µL of homogenization buffer.
- 3 To remove cell debris, samples were centrifuged at 2000 g for 5 min at 4 °C.
- 4 The supernatant was diluted 1:10 and 40 µL aliquoted into six sample wells of a 96-well plate.
- 5 At the time of assay, 160 µL of reduced cytochrome c (50 µM) was added to each sample well.
- 6 The plate was assayed in a Molecular Devices SpectraMax 384 Plus microplate reader, reading every 10 s at 550 nm for 5 min.
- 7 Oxidation of cytochrome c was indicated by a decline in the OD at 550 nm (extinction coefficient = 29.5).
- 8 Maximum slope of the optical density plot is calculated by the plate reader and is taken as V_{max} for the reaction.
- 9 A cytochrome c control was included and the protein concentration of each sample was determined using the BioRad DC Protein Assay method.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited